User Manual

Molded Fume Hood

Model : PMH-720, POH-720
Manual No. : B5011L001 Version : 2.2

⚠️ WARNING

Please read this user manual carefully before using the product. Users should operate the product according to the operating instructions contained in this user manual. Users should be aware of cautions and warnings for your safety. Failure to follow these guidelines may result in serious injury or damage to your property.

Thank you for purchasing JEIO TECH's products.

JEIO TECH Co., Ltd. is committed to offer the customer service including product sales, Technical service, and follow-up management. If you have any question about this user manual, please contact our Sales Department.

Quality Management System ISO 9001

ISO 9001 is a standard that defines the requirements for a Quality Management System (QMS). It ensures the suitability of implementation in order for suppliers to continuously provide products and services that meet the buyer's requirements. JEIO TECH Co., Ltd. obtained ISO 9001 certification from *Korean Foundation for Quality (KFQ)* since 1998 to provide better service. Under this certification system, we design, develop, produce and sell a wide range of products related to chemical engineering, environment and life. We also provide reliability testing equipment for electrical and electronic engineering. (Dryer, incubator, constant temperature and humidity chamber, constant temperature water tank, low temperature circulation tank, heat exchanger, stirrer, etc.)

You can download the above image by visiting the following JEIO TECH site.
<http://www.jeiotech.com>

Disclaimer

The material in this manual is for the users' information only. The products it describes are subject to change without prior notice, due to the manufacturer's continuous development program.

Copyright

No part of this manual may be reproduced or transmitted in any form or by any means, including photocopying, recording, or using information storage and retrieval systems, for any purpose other than the purchaser's own use, without the express written permission of JEIO TECH Co., Ltd.

©2020. All Rights Reserved. JEIO TECH Co., Ltd.

Third-party product names or brand names mentioned in this user manual may be the trademarks of their respective owners. JEIO TECH has no connection with those companies and we do not endorse or recommend them directly nor indirectly.

Table of Contents

1.0 Safety	1
1.1 How to Use the User Manual.....	1
1.2 Safety Signs for Safe Use of the Product.....	1
1.3 Disclaimer of Warranties	2
1.4 Warnings for safe use of the product.....	3
1.5 Caution Statement for Safe Use of Products.....	3
2.0 Functional Description	4
2.1 Product Introduction.....	4
2.2 Features	4
2.3 Structure	5
3.0 Installation.....	6
3.1 Checking the fume hood	6
3.2 Preparations before Installation	6
3.2.1 Minimum required space	6
3.2.2 Precautions for Selecting Fan Specifications.....	6
3.2.3 Environmental Conditions for hood Installation.....	7
3.2.4 Pre-checking points	7
3.3 Connection of filter type hazardous gas removal device (EM-33) (Optional).....	8
3.4 Bottom plate (Optional).....	9
4.0 Operation.....	10
4.1 Preparation of Work	10
4.2 How to Use the Fume Hood	10
4.3 End of Task	11
5.0 Maintenance	12
5.1 Maintenance and Management	12
5.2 Cleaning the product.....	12
5.2.1 Cleaning regularly	12
6.0 Troubleshooting	13
6.1 Abnormality handling.....	13
6.1.1 Abnormal behavior	13
7.0 Appendix	14
7.1 Technical Specifications.....	14
7.2 Accessories.....	14
7.3. Disposal of Products	15
7.4 Warranty.....	16
7.4.1 Service Warranty Period	16

7.4.2 Warranty Exceptions.....	16
7.4.3 Service Request Method.....	16
7.4.4 Returning procedure of the product	16
7.5 Service and Technical Support	17

1.0 Safety

1.1 How to Use the User Manual

This user manual contains instructions on how to use the product properly and safely. Please make sure to read and be aware of the safety warnings and precautions before use. Also, be sure to keep the user manual close to the product so that you can always use when needed. Also, please place the user manual in a prominent place on the product body so that new users can know how to use it safely and properly when transferring or leasing the product.

1.2 Safety Signs for Safe Use of the Product

- (1) We use safety signs to help control a hazard. They ensure safe use of the product and to prevent the potential accidents or damages.
- (2) The safety signs have been classified by the degree of hazards.
- (3) Please read carefully and understand the contents of the safety signs.

Safety Signs	Contents
 DANGER	Indicates a hazardous situation which, if not avoided, will result in death, fatal injury, or property damage
 WARNING	Indicates a hazardous situation which, if not avoided, could result in death, serious injury, or property damage
 CAUTION	Indicates a hazardous situation which, if not avoided, may result in minor or moderate injury
NOTICE	Notifies you careful actions that are not relevant to your injury

1.3 Disclaimer of Warranties

(1) We do not accept any warranty or liability for this product except for the specified in the company's assurance terms.

(2) This company is not responsible for any case where a user or third party has damaged directly or indirectly due to an unpredictable defect or natural disasters.

1.4 Warnings for safe use of the product

WARNING

Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

Please be fully aware of all the warning labels listed in this user manual.

Do not remove and damage the warning labels.

Do not move the hood while it is operating.

Do not use the hood as a storage for chemicals or devices.

Do not store flammable gas or flammable materials inside the fume hood.

Do not install the fume hood in environments that may cause flammable gas leaks.

Do not use the fume hood in environments where there is a possibility of explosion due to organic gas evaporation.

Do not expose the fume hood to direct sunlight.

Do not expose the fume hood to direct heat sources.

Do not disassemble or modify the fume hood at any time other than those described in this user manual.

Do not use damaged fume hoods.

1.5 Caution Statement for Safe Use of Products

CAUTION

Indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.

- When using the fume hood, work 6 inches (150mm) back from the front door.
- Do not insert the operator's head into the fume hood.
- Always keep the duct or ventilation system clean.
- Do not place chemicals or other devices on top of the hood while cleaning.
- Do not wash the hood with a high acid detergent, but please use a soft cloth.
- Avoid contact with insecticides or combustible gases.

2.0 Functional Description

2.1 Product Introduction

This product is Molded Fume Hood, which can safely release fine dusts and harmful gases generated when handling chemicals or working related to the product, to maintain a pleasant laboratory environment.


The application of the laboratory is as follows.

- PCR Operations
- Slide Staining Operations
- Drum decanting
- Robotics Enclosures
- Graphic arts preparation
- Powder weighing and dispensing
- Solvent bonding, welding or decanting
- A demonstration experiment showing


2.2 Features

1. The fume hood is made in a compact size, so it is easy to move, efficient space utilization, and easy to install and use.
2. Because the front door is made of transparent polycarbonate, it is not only easy to check the internal condition when working, but it is also firm.
3. The height of the front door can be adjusted depending on the user's needs with total of 5 levels at 60mm intervals.
4. The main body is made of polyethylene material, therefore it is not only firm but also durable for chemical resistant.
5. Flexible ducts with a diameter of 150mm can also be used using a fan in the existing laboratory.
6. You can remove harmful gases by connecting them to their filter-type Harmful Gas Removal (EM-33) or AH Series.
7. The inner corners of the fume hood are rounded for easy cleaning.
8. The front and the rear windows of the POH-720 model are transparent which is easy to observe the experimental processes on the both sides, and is suitable for use in demonstration at schools.

2.3 Structure

[The front view of PMH-720 and POH-720]



[The rear window of POH-720]

- (1) Body part
- (2) Door Handle
- (3) Duct Cover
- (4) Duct Hall
- (5) Tightening bolts for accessories
- (6) Door
- (7) Rear Window (POH-720 only)

3.0 Installation

3.1 Checking the fume hood

Please check the components listed in this User Manual after disassembling the package. If your order does not match the model name you ordered, or if you have not received any of the components below, please contact your JEIOTECH sales department or dealer.

Item	Figure	Quantity	Description
Main Body		1	-
User Manual		1	-

3.2 Preparations before Installation

3.2.1 Minimum required space

Keep at least 30 cm distance from other experimental products.

 WARNING

- Avoid heat sources or direct sunlight and place them on a flat floor.
-

3.2.2 Precautions for Selecting Fan Specifications

Do not use a fan with the inflow velocity between 0.3m and 1m/s.

3.2.3 Environmental Conditions for hood Installation

Install the Molded Fume Hood in a suitable space under the following conditions:

Avoid direct sunlight.

Keep the ambient temperature between 2 and 60°C.


Keep the relative humidity below 80%.

3.2.4 Pre-checking points

- Install the product in a flat place and do not cause any strong shock or vibration.
- Do not wash it with acidic detergent and please use a soft cloth.
- Do not place chemicals or other devices on top of the device while cleaning.
- Avoid contact with insecticides or combustible gases.
- Do not use the product in environments where there is a possibility of flammable gas leaks.
- Do not use the hood in environments where there is a possibility of explosion due to organic evaporation gas.
- Do not place explosive, flammable chemicals (alcohol, benzene, etc.) on the fume hood.

3.3 Connection of filter type hazardous gas removal device (EM-33) (Optional)

[Connecting the Filter-type Hazardous Gas Removal Device (EM-33) to the PMH-720/POH-720 Hood]

- (1) Combine EM-33 and PLINTH-PMH as indicated in the figure and tighten the four M4 bolts from the under direction.
 - (2) Place the plinth-equipped EM-33 model to the top of PMH-720 in place of the duct hole.
 - (3) Loosen the M6 plastic nut on the upper part of PMH-720 and tighten the M6 plastic bolt with PLINTH-PMH in the inner direction.
- ※ The POH-720 applies the same procedures.


3.4 Bottom plate (Optional)

You can choose polyethylene material for a bottom plate as an option that has higher chemical resistance when storing reagents.

When using a microscope or balance, a flat floor is a very important factor, so it is recommended that you purchase a bottom plate and use it.

Each bottom plate is equal in left and right side, so you can use on any side of the product.
Two sheets fit the product but just one sheet is also able to purchase and use on it.
(Refer to 7.2 Accessories)

It is separated into two sheets, so you can easily take only the contaminated plate out and clean it conveniently when there is a spillage. The corners are rounded for safer use.

[The appearance of bottom plate]

4.0 Operation

4.1 Preparation of Work


- (1) Check the connection between the fume hood, the duct and the ventilation system of the building.
- (2) Check the operation of the fan connected to the hood before use and run it.

4.2 How to Use the Fume Hood

You can adjust the height of the front door according to the following orders.

- (1) When lowering the door, use both hands to push down the handle on both ends of the door to the inner direction at the same time.
- (2) When raising a door, hold the door or handle and raise it.
- (3) Place the door handle in the desired position and make sure that it is correctly fixed.

In case of using the fan in 7.1, use the location at the bottom of the door, as shown below, within range of H1 to H2.

At the location above H2, the airflow velocity will be outside the reference range, which can make contaminants to leak to the operator and the laboratory.

4.3 End of Task

- (1) Stop the operation of the fan and the exhaust system.
- (2) Close the door downward as described in the procedure 4.2.

⚠ CAUTION

- Make sure that the front door handle is correctly fixed before operation.
 - Do not put flammable substances inside the fume hood.
 - When operating the fume hood, make a distance with 6 inches (150 mm) away from the front door.
 - Do not insert the operator's head inside the fume hood.
 - Do not remove the door of this fume hood.
 - Ensure that there are no combustion and explosive substances inside the hood.
 - Do not use this hood for storing chemicals or other devices.
 - Always keep the duct or building ventilation system clean.
-

5.0 Maintenance

5.1 Maintenance and Management

Category	Inspection Cycle		
	Everyday	Every week	Every year
Status of the duct connection with the product	•		
Surface damage and cleaning status of the product	•		
Cleaning status of door and duct	•		
Cleaning status of the product inside	•		
Check the building ventilation system		•	

5.2 Cleaning the product

5.2.1 Cleaning regularly

You must clean the product regularly to keep your product more efficiently and increase its life. Please clean the product once a week according to the procedure below.

Step 1: Remove chemicals or devices inside the product.

Step 2: Wipe lightly with a mild detergent and dry cloth.

WARNING

- Do not place chemicals or other devices on top of the product while you clean it.
- Do not use sulfuric acid, hydrochloric acid or organic solvents to clean the surface of the product.

CAUTION

- Do not put your head inside the product when cleaning it.
- Check if there is any foreign substance before use the product and keep them clean all the time.
- Do not allow contact with the body to high concentrations of nitric acid, sulfuric acid, sodium hydroxide, or corrosive solvents such as acetone, benzene, phenol, toluene, chloroform, crezole, acetic acid, chlorine, etc.

6.0 Troubleshooting

6.1 Abnormality handling

In case of abnormal use of the product, check the following information and follow the instructions. If you encounter problems that are not listed below, or if you can't resolve them by checking the following information, please contact your JEIOTECH sales department immediately to resolve the problem.

6.1.1 Abnormal behavior

Possible Problems	Cause types	Problem Solving
Front door handle is broken	<ul style="list-style-type: none">• Timeworn	<ul style="list-style-type: none">• Replace the hood parts.
Air discharge is poor	<ul style="list-style-type: none">• Poor connection between the product and duct• Abnormality of duct or building ventilation system• Inappropriate degree of door opening and closing• Fan failure	<ul style="list-style-type: none">• Connect duct with the product accurately.• Check the ventilation system of the duct or building.• Adjust the degree of opening and closing of the door within the appropriate range of use.• Replace the fan.
The door is damaged	<ul style="list-style-type: none">• Timeworn	<ul style="list-style-type: none">• Request for a technical support.




7.0 Appendix

7.1 Technical Specifications

Item / Model		Molded Fume hood / PMH-720	Molded Fume hood / POH-720
Dimension	Overall (W x D x H, mm)	720 x 472 x 535	
	Entrance (W x D, mm)	560 x 330	
	Rear window (W x D, mm)	-	604 x 380
	Duct (O.D, mm)	150	
Material	Body	Polyethylene	
	Door	Polycarbonate	
	Guide	Polypropylene	
Color	Body	White	
	Door / Rear window	Transparent	
	Guide	White	
Permissible air velocity (m/s)		0.3 ~ 1	
Weight (kg)		12	12
Cat. No.		AAAB5011	AAAB5021
Permissible ambient temperature (°C)		2 to 60	
Maximum relative humidity (%)		80	

※ The above specifications can be changed to improve quality and performance improvement.

7.2 Accessories

Cat. No.	Figures	Descriptions	Dimensions (W x D x H, mm)
AAAF1536		EM-33	348 X 498 X 343
AAAF15363		PLINTH-PMH	480 X 342 X 35
AAAB6411		Brackets for Molded Fume Hood	178 X 178 X 2
AAAF61501		Bottom plate	340 x 390 x 8

7.3. Disposal of Products

Before disposing a product or its component, please be aware of the followings:

1. This product is potentially hazard when exposed to biological, chemical, or radioactive materials, so it must be cleaned up to protect disposal personnel, waste collectors, and environment.
2. Please contact a person who is in charge of disposing of the products to check the electrical, electronic and material disposal standards. When you request externally, it is recommended that you contact to your company's vendor or a testing equipment manufacturer that can properly dispose the laboratory equipment and its components.

7.4 Warranty

7.4.1 Service Warranty Period

The warranty period is two years from the date of purchase, when a failure occurs while the device is in use. Within the service warranty period, you will be able to receive the repair service, according to the JEIOTECH warranty policy.

7.4.2 Warranty Exceptions

Paid repairs will be carried out even within the warranty period in the following cases:

- In case of failure due to natural disasters such as fire, flood, etc.
- In case of a failure due to dropping off or impacting on the product.
- Damage to the exterior caused by organic solvents such as excitement and benzene.
- Failure to comply with safety precautions and warnings within this user manual.
- In the event of a failure due to repair or modify the product by a person who is not delegated or not a JEIOTECH employee.

7.4.3 Service Request Method

If you experience any problems with your product, complete the service application form below and submit it to your head office or dealer for prompt and accurate repair.

- Purchase date
- Customer name / address / telephone number / e-mail
- Description of failure status

7.4.4 Returning procedure of the product

If you return your product for transportation-related damage and other reasons, you will receive your returning certificate number from your distributor. When returning your device, please complete the returning application from with the following information and send it to the headquarters or the dealer you purchased.

- Purchase date
- Customer name / address / telephone number / e-mail
- Delivery company's information
- Returning authentication number
- Reason for returning

7.5 Service and Technical Support

International Sales Head Office (Korea)

#1005, Byucksan Digital Valley 6-Cha, 481-4 Gasan-Dong, Geumcheon-Gu, Seoul 153-704, Korea

Tel: +82 2 2627 3816 E-mail: overseas@jeiotech.com
FAX: +82 2 3143 1824

The Americas (U.S.A. Branch)

1-A Gill St. Woburn, MA 01801, U.S.A.
Tel: +1 781 376 0700 E-mail: info@jeiotech.com
FAX: +1 781 376 0704

Europe (U.K. Branch)

Unit 3, Tower Industrial Park, Chalgrove, Oxfordshire, OX44 7XZ, United Kingdom
Tel: +44 1865 400321 E-mail: labcompanion@medlinescientific.com
FAX: +44 1865 400736


China (Shanghai Branch)

A-2113 Oriental International Plaza, 85 LouShanGuan Rd, Changning District, Shanghai, China 200336
Tel: +86-21-5108-9161 +86 21 3251 1086 E-mail: longjuncao@jeiotech.com
FAX: +86 21 3251 1083

South East Asia (Malaysia Branch)

No 57-59, Jalan Adenium 2G/6, Pusat Perniagaan Adenium, 48300 Bandar Bukit Beruntung, Selangor Darul Ehsan, Malaysia
Tel: +60 3 6028 5833 +60 3 6028 5825 E-mail: labcomp@streamyx.com
FAX: +60 3 6028 5822

- The contents of this manual are subject to change without notice to improve the performance of the product and to ensure accurate information delivery.
- Copyright of this manual is in JEIOTECH.

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©2020. All Rights Reserved. Jeio Tech Co., Ltd.