Nitro[®] Neuro[™] upright walker
Andador vertical Nitro[®] Neuro[™]
Déambulateur Nitro[®] Neuro[™]

item # / N° d'article / artículo #
102663RD

www.drivemedical.com

parts list

WEIGHT CAPACITY: 360 LBS (163 KG)

DO NOT use this product if your body weight exceeds the maximum weight capacity listed for this product.

- | | |
|----------------------------------|--------------------------------------|
| 1. Reflector | 11. Seat |
| 2. Hand Grip | 12. Seat Handle |
| 3. Brake Handle | 13. Storage Bag |
| 4. Padded Armrests | 14. Rear Caster |
| 5. Hand Grip Bar Adjustment Knob | 15. Brake Pad |
| 6. Accessory Mount | 16. Resistance Brake |
| 7. Height Adjustment Button | 17. Resistance Brake Adjustment Knob |
| 8. Sit to Stand Handles | |
| 9. Backrest | 18. Bag Mounts (2 locations) |
| 10. Backrest Adjustment Knob | |

Nitro Neuro Upright Walker includes:

- | | |
|--------------------------|---------------------------|
| 1. Removable Storage Bag | 3. Operating Instructions |
| 2. Backrest | 4. Accessory Mount |

unfolding the Nitro Neuro Upright Walker

To unfold the upright walker, position the Nitro Neuro on all four wheels and ensure all four wheels rest firmly on the ground. Hold each sit to stand handle and pull outward until the Nitro Neuro upright walker unlocks and opens fully (Figure 1). Gently push down on each side of the seat rails near the vertical frame to ensure the unit is fully opened. An audible 'click' confirms that the latches are secure and the seat is fully engaged in a locked position (Figure 2).

assembly instructions

1. Insert the 2 hand grip bars into the designated slots under the padded armrests. Rest your forearms naturally on the padded armrests so they sit flat and relaxed. (Figure 3)
2. Hand grip bar can be moved forward and back to a comfortable length. (Figure 4)
3. The hand grips should be within easy reach, allowing you to hold them without stretching or straining. (Figure 5)
4. The angle of the brake handles can be adjusted left and right to make your grip more comfortable.

assembly instructions continued

- Once you have found a comfortable position, secure the hand grip bar in place by tightening the hand grip bar adjustment knob underneath the tube clockwise. (Figure 6)

This upright walker features a reverse braking mechanism, meaning the brakes are always engaged by default to prevent unintended movement. To move the upright walker, the user must pull or squeeze the brake release handles.

Note: Pulling or engaging one or both brake release handles will unlock or release the brakes, allowing the upright walker to roll freely.

height adjustment

To adjust the handle height, press the height adjustment button and pull the hand grip upward at the same time (Figure 7). Release the button once the desired height is reached. An audible 'click' confirms the handle tube is secure in the locked position. Repeat this process until the desired height is reached for both handles. Refer to the numbering system (Figure 8) to ensure both handles are adjusted to the same height.

backrest adjustment

Turn backrest adjustment knob counterclockwise to loosen. Slide backrest to desired height then turn knob clockwise to tighten. (Figure 9)

DO NOT HANG ANYTHING ON THE BACKREST

walking and braking

The Nitro Neuro upright walker will not roll until you release the brakes.

1. Squeeze the left, right or both brake handles to release the brakes. (Figure 10)
2. Release brake handle(s) to stop. An audible click will indicate that the brakes have been activated upon releasing the brake handle(s).

If more or less resistance is needed to control walking speed, follow these steps to adjust the resistance brake.

1. Locate the resistance brake adjustment knob near the wheel – this controls how much friction is applied to the wheel when walking. (Figure 11)
2. To increase resistance (slower movement), turn the knob clockwise to add additional resistance to the wheel. Adjust both resistance brake knobs to same setting. (See Figure 12)
3. To decrease resistance (faster movement), turn the knob counterclockwise to reduce additional resistance to the wheel. Adjust both resistance brake knobs to same setting. (See Figure 12)

removing / attaching the storage bag

To attach the storage bag to the upright walker, affix the bag mounting clips (located on either side of the storage bag) to the desired mounting points on the upright walker (Figures 13 and 14). Easily remove the storage bag and attach the shoulder strap using the clips located on either side of the storage bag to take it with you on the go. (Figure 15)

**WEIGHT CAPACITY OF THE BAG IS 11 LB.
DO NOT OVERLOAD AS IT WILL AFFECT PERFORMANCE.**

sitting

While facing the upright walker, position the open rollator on a flat, even, solid surface on all four wheels. The brakes are always engaged unless you pull the handle brake release. The upright walker must remain upright and must not move. Utilize the sit to stand handles for support, to turn completely around. With both hands on the handles slowly sit, ensuring you sit straight up and do not lean backwards. (Figure 16).

DO NOT MOVE UPRIGHT WALKER WHILE SEATED.

folding / transport / storage

To fold, pull the seat upward by the seat handle (Figure 17). The upright walker can now be folded completely and secured by pushing the sit to stand handles together. An audible 'click' confirms that the under-seat locking mechanism is engaged and the Nitro Neuro is locked in the folded position (Figure 18).

To transport the rollator, place one hand on the outside of each frame. Lift carefully using proper lifting technique: keep your back straight, bend at the knees, and use your legs to lift.

The rollator weighs approximately 25 lb. Use caution when lifting.

cleaning

Clean regularly by hand with a damp cloth and mild detergent only.

Do not use aggressive detergents such as chlorine, bleach or abrasive cleaners.

maintenance

Drive Medical recommends periodic visual inspection of this product by the user to ensure that all parts and hardware are secure and all components are in good working order and not worn, torn, frayed, or loose. There should be no obstructions or impediments to normal and safe operation. If any part or mechanism appears to be unsafe or compromised, do not use the product. If you experience any problems or issues, please contact the provider that sold you the product or visit www.drivemedical.com/help for more information.

Periodically inspect the moving components of the Nitro Neuro Upright Walker for signs of wear or damage. On a daily basis, test the braking system to ensure it is functioning properly.

If you experience any issues with wheel resistance, braking performance, or the resistance brake adjustment system, please contact your authorized Drive dealer or Drive customer service for assistance.

Important: Do not pull on the brake cables. Pulling on the cables may cause them to kink, stretch, or become damaged, which could prevent the braking system from operating correctly. Damaged cables should be replaced immediately. If the cabling becomes damaged, discontinue use and contact your authorized Drive dealer or Drive customer service for assistance.

safety instructions

- This is a walking aid only and is not to be used as a transportation device.
- Use only on even and solid surfaces.
- To prevent unintended movement, always maintain control of the upright walker
- The Nitro Neuro upright walker will not roll until you release the brakes.
- Pulling or engaging either brake handle(s) will disengage the brake, allowing the upright walker to roll freely. Ensure the user is prepared and stable before releasing the brake.
- Maximum user weight: 360 lb.
- Regularly check the screws and fasteners on the Nitro Neuro upright walker and retighten if necessary.
- Use bag only for transport of objects and do not overload. Maximum load: 11 lb.
- Do not expose the product to temperatures above 104°F for an extended period as this can cause the Nitro Neuro upright walker parts to overheat and burn the user.

- Please make sure the upright walker is fully opened, seat is down and the folding bar under the seat is locked in the down position before using or sitting in the rollator.
- Be sure both hand grips are at equal heights. Test upright walker for stability before applying full weight to sit to stand handles or armrest and handles.
- Replace casters when significantly worn or damaged.
- The backrest is designed to provide postural support while seated. It is not intended to support the user while standing, or be used to hang any items.
- **DO NOT MOVE UPRIGHT WALKER WHILE SEATED.**
- Consult your physician or healthcare professional before using this product to ensure it is appropriate for your individual needs.

lifetime limited warranty

Your Drive® brand product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence or normal wear and tear. The warranty does not extend to non-durable components such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive® device or this warranty, please contact an authorized Drive® dealer.

WARNING: This product can expose you to chemicals, including di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

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© 2026 Medical Depot, Inc. All rights reserved.
Drive and Nitro are registered trademarks of Medical Depot, Inc.
Sprint is a trademark of Medical Depot, Inc.
© 2026 Medical Depot, Inc. Tous droits réservés.
Drive et Nitro sont des marques déposées de Medical Depot, Inc.
Sprint est une marque de commerce de Medical Depot, Inc.