

Product Name:	Multi-layer Cell Culture Flask		
Product Application:	NEST cell culture flasks are an ideal solution for enhanced cell attachment and growth of a variety of primary cells and transformed cells in serum free or serum-containing cultures.		
Product Materials:	PS (bottle) + PP (cap) Meets USP, Class VI standards		
product Specifications:	Temperature range: stored at room temperature. Shelf life: 3 Years after date of production (ensure package is in good)		
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137		
Pyrogens:	Non-Pyrogenic		
RNase/DNase Testing:	DNase/RNase free		
BSE/TSE:	These products are deemed animal free.		
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Flatness testing: qualified Hydrophilic testing: qualified Packaging inspection: qualified		
Features:	1. High clarity medical grade polystyrene material 2. Two types of cap style: vent caps & plug seal cap 3. Vent caps with 0.22um hydrophobic filters for gas exchange without contamination 4. Culture surface TC treated, good for cell attachment and growth. 5. Growth area: 870 cm² 6. Individually Packaged in sterile bag.		
Guidelines For Use:	Culture of cells: 1. Mix cell suspension with medium: Prepare cell suspension of required concentration in a container. Recommended volume is about 30- 50mL per layer. 2. Add the mixed liquid into the Multi-layer Flask slowly with serological pipettes. To avoid foam and bubbles, allow liquid stream to flow along the slope of the Multi-layer		