

ABS



American BioTech Supply

PH-ABT-HC-UCFS-0504G

Premier Pharmacy/Vaccine Refrigerator



PH-ABT-HC-UCFS-0504G Pharmacy Refrigerator

Premier Undercounter Refrigerator (Free Standing)

5.2 Cu. Ft. Capacity

One glass swing door, right hinged, not reversible

Microprocessor temperature controller with battery back-up (batteries included)

Digital temperature display

Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)

Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration

Audible and visual high and low temperature alarms with battery back-up

Remote alarm contacts

Three adjustable shelves

Bottom footed shelf to allow air circulation

Back wall evaporator cover

Forced draft circulation

Keyed door lock

Magnetic door gasket for positive seal

Leveling legs

LED interior light

Switch activated light

Probe access port (3/8")

High density urethane foam cabinet insulation

HFC-free refrigerant (R600a), foam insulation and packaging

Pharmacy refrigerator/freezer toolkit and temperature logs

Two year parts and labor warranty

Exterior dimensions: 23 3/4"W x 24 1/4" D x 32 1/8" H

Note: This unit must have 4" clearance on sides and back for adequate ventilation

Cycle defrost

115V, 60 Hz, 1.3 Amps, 1/6 HP

UL/C-UL listed