

Advanced, multi user Bluetooth® blood pressure monitor with enhanced features that provide trusted accuracy and reliable measurements every time.

Product Benefits


- Experience accurate readings with advanced cuff fit detection, avoiding errors from incorrect cuff placement
- Ensure precise measurements with integrated movement detection, minimizing the risk of inaccurate results due to movement
- Automatically records and tracks 60 measurements with the highly rated Heart Track app to enable trending and sharing for better health management
- Displays the average of the last three measurements providing more consistent blood pressure data

Comfort and Convenience

- Up to 4 users provides convenience, cost savings, and enhanced data management capabilities for households with multiple individuals needing to monitor their blood pressure regularly
- Advanced inflation and deflation technology provides faster and more comfortable readings

The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value



Advanced Averaging

60 Measurement Memory Per User

Movement Error

Pressure Rating Indicator

Battery Life Indicator

Multi User



Irregular Heartbeat Detection (IHB)

Cuff Fit Error

One Button Operation

Unlock the Potential
Download the FREE
A&D Heart Track
Companion App



Technical Specifications

Measurement Method	Oscillometric
Memory	60 measurements per user
Cuff Size	8.6-16.5" (22-42 cm)
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats/min
Measurement Accuracy	Pressure: ±3 mmHg, Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	5.5 x 4.6 x 2.5 in (14 x 11.8 x 6.4 cm)
Monitor Weight	10.4 oz (295g) without batteries
Display Size (WxH)	1.7 x 1.9" (4.3 x 4.9cm)
Clinical Test	ISO81060-2: 2020
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Case Cube	UPC / GTIN
Individual Unit Box	5.1 x 4.1 x 7.7 in (13 x 10.5 x 19.5 cm)	1.7 lb (750 g)	0.09 ft ³ (0.003 m ³)	0-93764-60450-8
Case (10 units)	11 x 21.7 x 8.7 in (28 x 55 x 22 cm)	18.1 lb (8.2 kg)	1.2 ft ³ (0.03 m ³)	1-00-93764-60450-5

Optional Accessories

Name	Order Code
AC Power Adapter	TB-233

Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 6.3-9.4" (16-24 cm)	UA-289A
Medium Cuff 9.0-14.6" (23-37 cm)	UA-290A
Large Cuff 12.2-17.7" (31-45 cm)	UA-291A
Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-420A
Cuffs are not made with natural rubber latex	

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff 8.6-16.5" (22-42 cm)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 4 AA batteries