

ADA305AFIF



32" x 20.25" x 23" (H x W x D)

Built-in undercounter -25°C ADA compliant all-freezer with panel-ready door, factory installed lock, digital controls, interior baskets, hospital cord with 'green dot' plug, access port, and manual defrost operation

Highlights:

Flexible design allows built-in or freestanding use in 21" wide spaces

Easy temperature management with a digital thermostat

Basket shelving offers convenient storage of a variety of items

Product Features:

| | |
|---|--|
| 20" Wide Built-In Healthcare All-Freezer, ADA-Compliant | Slim-fitting all-freezer for general healthcare applications at an ADA-compliant height |
| ADA-Compliant Design | 32" height allows easy installation under lower ADA compliant counters |
| -25°C Capable | Low temperature operation ideal for vaccines, pharmaceuticals, and other sensitive materials |
| Built-In Capable | Front-breathing design lets you make the best use of space by installing your appliance under the counter |
| Panel-Ready Door | Integrated stainless steel door frame accepts full overlay panels, allowing you to customize the refrigerator door or blend it into your cabinetry design; panel size:19 13/16" W x 27 1/8" H; customers provide their own handles |
| Factory-Installed Lock | Located towards the top of the door for easy access and to provide added security of unit (2 keys included) |
| Access Port | Unit comes with a factory installed probe hole for users to connect their own monitoring equipment |
| Adjustable Digital Thermostat | Electronic controls located in the front kickplate for easy & intuitive temperature management, with a readout in Fahrenheit or Celsius |

ADA305AFIF Specifications:

| Overview | |
|----------------------------|--------------------|
| Nominal Height | 32.0" (81 cm) |
| Height of Cabinet | 31.75" (81 cm) |
| Height to Hinge Cap | 32.0" (81 cm) |
| Nominal Width | 21.0" (53 cm) |
| Width | 20.25" (51 cm) |
| Nominal Depth | 23.0" (58 cm) |
| Depth | 23.0" (58 cm) |
| Depth with door at 90° | 40.25" (102 cm) |
| Capacity | 2.47 cu.ft. (70 L) |
| Defrost Type | Manual |
| Door | Panel-Ready |
| Cabinet | White |
| US Electrical Safety | ETL |
| Canadian Electrical Safety | ETL-C |
| Amps | 1.1 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 75.0 lbs. (34 kg) |
| Parts & Labor Warranty | 1 Year |
| Compressor Warranty | 5 Years |
| UPC | 761101151298 |
| Freezer | |
| Door Swing | RHD |
| Reversible | Factory Reversible |
| Interior Height | 23.75" (60 cm) |
| Interior Width | 14.5" (37 cm) |
| Interior Depth | 17.0" (43 cm) |
| Shelf Type | Drawers |
| Shelf Qty | 3 |
| Thermostat Type | Digital |
| Fan Type | Exterior |
| Refrigerant Type | R600a |
| Refrigerant Amount | 1.34 oz. |
| High Side PSI | 180.0 |
| Low Side PSI | 80.0 |
| Level Legs | 2 |
| Interior Drawers | Yes |
| Drawers Qty | 3 |

| | |
|------------------------|---------------|
| Compressor Step Height | 5.75" (15 cm) |
| Compressor Step Width | 14.5" (37 cm) |
| Compressor Step Depth | 6.0" (15 cm) |
| Temperature Range | -25° to -15°C |

