

AGP34RFADA



32.75" x 19" x 19.75" (H x W x D)

General purpose ADA height 2-door refrigerator-freezer with digital refrigerator control, internal fan, and alarm/thermometer

Highlights:

Combination refrigerator-freezer for general purpose storage in facilities

Digital controls allow you easily manage the refrigerator temperature between 2 and 8°C

Alarm/thermometer included to help monitor the interior temperature

Product Features:

Refrigerator-freezer combination	2 to 8°C refrigerator section for basic storage applications, with an upper freezer compartment ideal for ice packs
Lower height for ADA applications	Height to cabinet is 32 5/8", with an additional 1/8" required for the hinge cap
Digital thermostat (refrigerator)	Easy operation with a factory-installed digital thermostat that allows you to manage the refrigerator temperature between 2 and 8°F (NOTE: the freezer temperature will self-adjust based on the refrigerator setting)
Alarm/thermometer included	Temperature monitoring device (TMD) displays the current and min/max temperature of the refrigerator interior, with an audible alarm to alert if the unit goes out of range
Buffered temperature probe	Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of refrigerator air temperature
External temperature readout	Thermometer is designed to be placed on the door, allowing you to view the interior temperature without opening the door and exposing stored contents to warmer air
Fan-forced cooling	Refrigerator section includes a factory-installed fan to help ensure a more evenly cooled interior
Automatic defrost refrigerator	Lower refrigerated section does not require manual defrosting

Manual defrost freezer Static cooling system in the upper freezer compartment

AGP34RFADA Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	32.63" (83 cm)
Height to Hinge Cap	32.75" (83 cm)
Nominal Width	19.0" (48 cm)
Width	19.0" (48 cm)
Nominal Depth	20.0" (51 cm)
Depth	19.75" (50 cm)
Depth with door at 90°	36.75" (93 cm)
Capacity	3.2 cu.ft. (91 L)
Defrost Type	Cycle
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	55.0 lbs. (25 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101086026
Refrigerator-Freezer	
Door Swing	RHD
Reversible	Yes
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	18.5" (47 cm)
Refrigerator - Interior Width	15.75" (40 cm)
Refrigerator - Interior Depth	15.0" (38 cm)
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	2
Refrigerator - Full Door Shelves	2
Refrigerator - Temperature Range	2 to 8°C
Freezer - Interior Height	7.5" (19 cm)
Freezer - Interior Width	14.0" (36 cm)
Freezer - Interior Depth	14.0" (36 cm)

Freezer - Full Door Shelves	1
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	1.16 oz.
High Side PSI	250.0
Low Side PSI	88.0
Level Legs	2
Compressor Step Height	9.0" (23 cm)
Compressor Step Width	15.75" (40 cm)
Compressor Step Depth	5.0" (13 cm)

