

AWC6PBLHD



35.5" x 23.5" x 23.88" (H x W x D)

Performance Series hospital blanket warmer with a self-closing glass door, stainless steel interior, digital thermostat, casters, and lock

Highlights:

Factory installed lock allows you to restrict access to staff

Full length glass door keeps complete interior in view for easy inventory check

Antimicrobial handle to help reduce germ spread

Product Features:

| | |
|--|---|
| Performance Series Blanket Warmer | Designed to keep blankets warm for patient care in hospitals, doctors' offices, physical therapy centers, and similar settings |
| Pre-set to 130°F | Standard temperature for keeping blankets warm for patient care (NOTE: parameters are locked at a 130°F maximum, but can be adjusted up to 160°F. Contact the Accucold technical support team for instructions) |
| Microprocessor control panel | Digital controls with a large LED readout in Celsius or Fahrenheit for intuitive temperature management |
| Antimicrobial silver ion-coated handle | Powder-coated handle with naturally occurring silver ions to help reduce the spread of germs |
| Double pane tempered glass door | Heat-safe glass door offers a full view of stored contents |
| Self-closing door | Door will automatically close when left ajar, preventing items from accidentally being left exposed to ambient temperatures |
| Stainless steel interior | Interior walls are constructed from polished 304 grade stainless steel to ensure lasting durability and easy cleaning |
| Adjustable stainless steel shelving | Heavy duty capacity (up to 44 lbs.) with a perforated design to help ensure proper heat distribution |

Hospital grade casters | 7" high rubberized casters (2 locking) ensure easy mobility as needed

AWC6PBLHD Specifications:

| Overview | |
|----------------------------|--|
| Height of Cabinet | 35.5" (90 cm) |
| Width | 23.5" (60 cm) |
| Depth | 23.88" (61 cm) |
| Depth with Handle | 25.38" (64 cm) |
| Depth with door at 90° | 45.63" (116 cm) |
| Capacity | 5.22 cu.ft. (148 L) |
| Door | Glass |
| Cabinet | White |
| US Electrical Safety | ETL |
| Canadian Electrical Safety | ETL-C |
| Amps | 2.0 |
| Voltage/Frequency | 115 V AC/60 Hz |
| Weight | 114.0 lbs. (52 kg) |
| Parts & Labor Warranty | 1 Year |
| UPC | 761101155166 |
| Warming Cabinet | |
| Door Swing | LHD |
| Reversible | No |
| Adjustable Shelves | Yes |
| Shelf Type | Steel |
| Shelf Qty | 4 |
| Thermostat Type | Digital |
| Interior Light | Yes |
| Interior Height | 23.5" (60 cm) |
| Interior Width | 19.75" (50 cm) |
| Interior Depth | 18.0" (46 cm) |
| Wheels | 4 |
| Temperature Range | 90 to 130°F; To adjust up to 160°F, contact our technical support team |

