

BCW100



27.25" x 16.75" x 15.25" (H x W x D)

Glass door warmer to hold 54 packs of disposable bathing wipes

Highlights:

Stores up to 54 8-packs of disposable bathing cloths

Designed to operate at 125°F

Versatile design can be used freestanding, tabletop, or wall-mounted

Product Features:

Bathing Wipe Warming Cabinet	Designed to keep disposable wipes warmed at 125°F
Full capacity	Each pocket is designed to hold three sealed 8-piece packets of bathing cloths
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents, allowing staff to check inventory without opening the door
Reversible door	User-reversible door swing for added flexibility to best accommodate your setting
Quiet operation	No moving parts to ensure a quiet performance
Overheat protection	Helps to ensure interior temperature stays steady at 125°F
Stainless steel pocket racks	Convenient separation of bathing wipe packets with steel walls that are easy to clean
Leveling legs	Adjustable legs allow you to level the unit as needed
Flush back	Space saving design with easier cleanability

BCW100 Specifications:

Overview	
Height of Cabinet	27.25" (69 cm)
Width	16.75" (43 cm)
Depth	15.25" (39 cm)
Depth with Handle	15.88" (40 cm)
Depth with door at 90°	31.0" (79 cm)
Capacity	0.65 cu.ft. (18 L)
Door	Glass
Cabinet	Black
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	3.04
Voltage/Frequency	115 V AC/60 Hz
Weight	44.0 lbs. (20 kg)
Parts & Labor Warranty	1 Year
UPC	761101161204
Warming Cabinet	
Door Swing	RHD
Reversible	Yes
Shelf Qty	18
Interior Height	22.0" (56 cm)
Interior Width	14.25" (36 cm)
Interior Depth	11.5" (29 cm)
Level Legs	4
Temperature Range	125°F

