

FF28LWHPLUS2



24.88" x 18.5" x 17.75" (H x W x D)

Compact all-refrigerator with automatic defrost, digital thermostat, internal fan, NIST calibrated thermometer, and lock; replaces FF28LWHPLUS

Highlights:

- NIST calibrated thermometer displays the current and high/low temperature to the nearest tenth of a degree
- Factory installed lock provides security you can count on
- Internal fan-forced cooling for more temperature stability

Product Features:

Automatic defrost	Reduced maintenance with auto defrost system
Factory installed lock	Keyed lock for a secure interior
NIST calibrated temperature display	Thermometer provides continuous display of the current and high/low temperature to the nearest tenth of a degree, with a readout switchable from Celsius to Fahrenheit at the touch of a button
Buffered temperature probe	Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of refrigerator air temperature
Digital thermostat	Electronic controls ensure easy and intuitive temperature management
Internal fan	Fan-forced cooling ensures the unit is evenly cooled from top to bottom
Flat door liner	Interior design conforms to CDC guidelines to ensure stable conditions for sensitive materials
Wire shelving	Helps to ensure improved air circulation, with a full-width fixed wire shelf and a removable half shelf included
Compact dimensions	Ideal size for use on top of counters

FF28LWHPLUS2 Specifications:

Overview	
Nominal Height	25.0" (64 cm)
Height of Cabinet	24.5" (62 cm)
Height to Hinge Cap	24.88" (63 cm)
Nominal Width	19.0" (48 cm)
Width	18.5" (47 cm)
Nominal Depth	18.0" (46 cm)
Depth	17.75" (45 cm)
Depth with door at 90°	34.5" (88 cm)
Capacity	2.4 cu.ft. (68 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	0.7
Voltage/Frequency	115 V AC/60 Hz
Weight	36.0 lbs. (16 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101057514
Refrigerator Features	
Door Swing	RHD
Reversible	Factory Reversible
Shelf Type	Wire
Shelf Qty	2
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.7oz.
High Side PSI	120.0
Low Side PSI	70.0
Level Legs Qty	2
Temperature Range	35°F to 50°F; for a tighter range, contact our product specialists

Dimensions

Interior Height	21.5" (55 cm)
Interior Width	16.0" (41 cm)
Interior Depth	13.5" (34 cm)
Compressor Step Height	7.75" (20 cm)
Compressor Step Width	16.0" (41 cm)
Compressor Step Depth	5.0" (13 cm)

