

# FF6LWPRO



33.5" x 23.63" x 23.5" (H x W x D)

24" wide auto defrost all-refrigerator with lock, digital thermostat, internal fan, door storage, and access port for user-provided monitoring equipment

**Highlights:**

Ideal 24" footprint offers large capacity in a slim fit

Factory installed access port allows users to install their own monitoring equipment

Advanced precision with digital temperature controls

**Product Features:**

|                                           |                                                                                             |
|-------------------------------------------|---------------------------------------------------------------------------------------------|
| Automatic defrost                         | Reduced maintenance with auto defrost system                                                |
| Digital thermostat                        | Electronic controls ensure easy and intuitive temperature management                        |
| Internal fan                              | Fan-forced cooling ensures the unit is evenly cooled from top to bottom                     |
| Access port                               | Includes a factory installed probe hole for users to install their own monitoring equipment |
| Hospital grade cord with 'green dot' plug | Secure three-pronged cord for added unit safety                                             |
| Factory installed lock                    | Top-mounted keyed lock for convenient security                                              |
| Adjustable shelves                        | Rearrange refrigerator space to accommodate all sizes or remove shelves for simple cleanup  |
| Door storage                              | Store items right on the door for easy access                                               |
| Interior light                            | Automatically illuminates when you open the door                                            |

## FF6LWPRO Specifications:

| Overview                   |                       |
|----------------------------|-----------------------|
| Nominal Height             | 34.0" (86 cm)         |
| Height of Cabinet          | 33.5" (85 cm)         |
| Nominal Width              | 24.0" (61 cm)         |
| Width                      | 23.63" (60 cm)        |
| Nominal Depth              | 24.0" (61 cm)         |
| Depth                      | 23.5" (60 cm)         |
| Depth with Handle          | 25.13" (64 cm)        |
| Depth with door at 90°     | 46.13" (117 cm)       |
| Capacity                   | 5.5 cu.ft. (156 L)    |
| Defrost Type               | Automatic             |
| Door                       | White                 |
| Cabinet                    | White                 |
| US Electrical Safety       | ETL                   |
| Canadian Electrical Safety | ETL-C                 |
| Amps                       | 1.4                   |
| Voltage/Frequency          | 115 V AC/60 Hz        |
| Weight                     | 100.0 lbs. (45 kg)    |
| Parts & Labor Warranty     | 1 Year                |
| Compressor Warranty        | 5 Years               |
| UPC                        | 761101070667          |
| Refrigerator Features      |                       |
| Door Swing                 | RHD                   |
| Reversible                 | Yes                   |
| Shelf Type                 | Wire                  |
| Shelf Qty                  | 3                     |
| Full Door Shelf Qty        | 3                     |
| Adjustable Shelves         | Yes                   |
| Thermostat Type            | Digital               |
| Fan Type                   | Interior              |
| Refrigerant Type           | R600a                 |
| Refrigerant Amount         | 0.64oz.               |
| High Side PSI              | 180.0                 |
| Low Side PSI               | 80.0                  |
| Level Legs Qty             | 2                     |
| Interior Light             | Yes                   |
| Temperature Range          | 2 to 8°C / 36 to 46°F |
| Dimensions                 |                       |

|                        |                |
|------------------------|----------------|
| Interior Height        | 27.75" (70 cm) |
| Interior Width         | 20.0" (51 cm)  |
| Interior Depth         | 17.5" (44 cm)  |
| Compressor Step Height | 7.0" (18 cm)   |
| Compressor Step Width  | 20.0" (51 cm)  |
| Compressor Step Depth  | 6.0" (15 cm)   |

