

FF7LWGPCADALHD



32.25" x 23.63" x 23.5" (H x W x D)

Commercially listed ADA height all-refrigerator for general purpose healthcare use with front lock

Highlights:

32" height for recessed use in settings complying with ADA guidelines for lower counter height

Commercially approved for use in foodservice establishments

Factory-installed lock provides security you can count on

Product Features:

24" Wide Healthcare Refrigerator	Designed for general purpose storage in institutional facilities
ADA compliant	Sized at 32" high for use in settings complying with ADA guidelines for lower counter height
Commercially approved	ETL-S listed to NSF-7 commercial standards
Fully finished white cabinet	Allows the refrigerator to be used freestanding
Front-mounted lock	Keyed lock for a secure interior
Automatic defrost operation	Reduced user maintenance with automatic defrost
Internal fan	Internal fan ensures the interior is evenly cooled from top to bottom
Mechanical thermostat	Dial thermostat located inside for easy temperature management
One piece seamless interior liner	Enjoy easy clean-up with a one piece liner that won't hold a mess in hidden crevices

FF7LWGPCADALHD Specifications:

Overview	
Height of Cabinet	32.0" (81 cm)
Height to Hinge Cap	32.25" (82 cm)
Width	23.63" (60 cm)
Depth	23.5" (60 cm)
Depth with Handle	25.13" (64 cm)
Depth with door at 90°	46.13" (117 cm)
Capacity	5.5 cu.ft. (156 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Amps	0.8
Voltage/Frequency	115 V AC/60 Hz
Weight	100.0 lbs. (45 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101161594
Refrigerator Features	
Door Swing	LHD
Shelf Type	Wire
Shelf Qty	3
Adjustable Shelves	Yes
Thermostat Type	Dial
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.64oz.
Level Legs Qty	2
Interior Light	Yes
Temperature Range	1 to 10°C
Dimensions	
Interior Height	27.75" (70 cm)
Interior Width	20.0" (51 cm)
Interior Depth	17.5" (44 cm)
Compressor Step Height	7.0" (18 cm)
Compressor Step Width	20.0" (51 cm)

Compressor Step Depth

