

FF7LWPRO



33.5" x 23.63" x 23.5" (H x W x D)

24" wide auto defrost commercial all-refrigerator with lock, digital thermostat, internal fan, and access port for user-provided monitoring equipment

Highlights:

Ideal 24" footprint offers large capacity in a slim fit

Factory installed access port allows users to install their own monitoring equipment

Commercially listed to NSF-7 standards

Product Features:

Commercially approved	ETL-S listed to NSF standards for commercial use
Digital thermostat	Electronic controls ensure easy and intuitive temperature management
Internal fan	Fan-forced cooling ensures the unit is evenly cooled from top to bottom
Access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Hospital grade cord with 'green dot' plug	Secure three-pronged cord for added unit safety
Factory installed lock	Top-mounted keyed lock for convenient security
Automatic defrost	Reduced maintenance with auto defrost system
Adjustable shelves	Rearrange refrigerator space to accommodate all sizes or remove shelves for simple cleanup
Flat door liner	Easy to clean, with more depth for storage

FF7LWPRO Specifications:

Overview	
Nominal Height	34.0" (86 cm)
Height of Cabinet	33.5" (85 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.5" (60 cm)
Depth with Handle	25.13" (64 cm)
Depth with door at 90°	46.13" (117 cm)
Capacity	5.5 cu.ft. (156 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Amps	1.4
Voltage/Frequency	115 V AC/60 Hz
Weight	100.0 lbs. (45 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101072111
Refrigerator Features	
Door Swing	RHD
Reversible	Yes
Shelf Type	Wire
Shelf Qty	3
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.64oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs Qty	2
Interior Light	Yes
Temperature Range	2 to 8°C / 36 to 46°F
Dimensions	

Interior Height	27.75" (70 cm)
Interior Width	20.0" (51 cm)
Interior Depth	17.5" (44 cm)
Compressor Step Height	7.0" (18 cm)
Compressor Step Width	20.0" (51 cm)
Compressor Step Depth	6.0" (15 cm)

