

FS407LW



34" x 19.25" x 21.38" (H x W x D)

20" wide counter height manual defrost all-freezer for general purpose healthcare use, -20°C capable with a lock and white exterior

Highlights:

-20° C capable with interior slide-out drawers for easy storage

Manual defrost operation for improved stability in healthcare facilities

Factory-installed lock provides security you can count on

Product Features:

General purpose healthcare design	Low temperature manual defrost freezer for use in healthcare, life science, and other institutional facilities
Slim 20" width	2.8 capacity inside a slim footprint
-20°C capable	Low temperature capability ideal for certain sensitive items
Factory installed lock	Located towards the top of the door for easy access and to provide added security of unit (2 keys included)
Manual defrost	Efficient performance with more reliable temperature stability
Basket shelves	Four slide-out baskets offer convenient storage, with plastic fronts to help protect contents from warmer air when the door is opened
Hospital grade cord with 'green dot' plug	Secure three-pronged cord for added unit safety
Adjustable thermostat	Rear mounted dial control to prevent accidental adjustments
Flat door liner	Easy to clean, with more depth for storage

FS407LW Specifications:

Overview	
Nominal Height	34.0" (86 cm)
Height of Cabinet	33.5" (85 cm)
Height to Hinge Cap	34.0" (86 cm)
Nominal Width	20.0" (51 cm)
Width	19.25" (49 cm)
Nominal Depth	22.0" (56 cm)
Depth	21.38" (54 cm)
Depth with Handle	23.0" (58 cm)
Depth with door at 90°	39.5" (100 cm)
Capacity	2.8 cu.ft. (79 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.6
Voltage/Frequency	115 V AC/60 Hz
Weight	57.0 lbs. (26 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101083551
Freezer	
Door Swing	RHD
Reversible	Yes
Interior Height	26.38" (67 cm)
Interior Width	14.0" (36 cm)
Interior Depth	16.75" (43 cm)
Thermostat Type	Dial
Refrigerant Type	R600a
Refrigerant Amount	1.23 oz.
High Side PSI	125.0
Low Side PSI	68.0
Level Legs	2
Interior Drawers	Yes
Drawers Qty	4
Temperature Range	-10 to -20°C

