

VLT650ADA



32" x 23.63" x 23.5" (H x W x D)

Built-in undercounter ADA compliant laboratory freezer capable of -35° C (-31° F) operation

Highlights:

Commercially listed to NSF-7 standards

Capable of operating at -35° C

Precise temperature control with a digital thermostat

Product Features:

ADA compliant	32" height allows installation under lower ADA compliant counters
Commercially approved	Meets NSF-7 commercial standards
-35° C (-31° F) capable	Low temperature operation ideal for vaccines, pharmaceuticals, and other sensitive materials
Built-in capable	Front-breathing system allows the unit to be used built-in
Factory installed lock	Keyed lock offers convenient security
Three slide-out drawers	Conveniently organize your items while keeping the temperature stable when the door is opened
Digital thermostat	Electronic controls for precise temperature management
Alarm with temperature display	Audible alarm sounds to warn of temperature variations with a digital readout on the thermometer to provide the interior temperature
Manual defrost	Static cooling system with manual defrost operation

VLT650ADA Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	32.0" (81 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.5" (60 cm)
Capacity	3.5 cu.ft. (99 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S
Voltage/Frequency	115 V AC/60 Hz
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101055077
Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Thermostat Type	Digital
Refrigerant Type	R134a
Refrigerant Amount	2.76 oz.
Interior Drawers	Yes
Drawers Qty	3

