

VT65MLBIVAC



33.5" x 23.63" x 23.63" (H x W x D)

Built-in undercounter vaccine storage all-freezer with digital thermostat, temperature alarm, hospital grade cord, self-closing door, and front lock

Highlights:

Medical grade design conforms with CDC guidelines for vaccine storage

Capable of operating at -25° C for long-term preservation of temperature sensitive materials

Interior drawers ensure stored items are not exposed to room temperature when the door is open

Product Features:

Slim undercounter dimensions	24" footprint is ideal for use under the counter
Built-in capable	Front-breathing design allows built-in installation under counters
-25° C/-13° F capable	Low temperature operation ideal for vaccines, pharmaceuticals, and other sensitive materials
Digital thermostat	Electronic controls for more precise temperature management, with temperature readout and controls externally located for added convenience
High/low temperature alarm	Audible alarm sounds if the interior temperature rises or falls out of set range
NIST calibrated temperature display	Alarm includes NIST calibrated display of the current interior and high/low temperature to the nearest tenth of a degree
Self-closing door	Door closes automatically when opened at angles less than 90°
3/8" access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Hospital grade cord and plug	Secure three-pronged cord for added unit safety

VT65MLBIVAC Specifications:

Overview	
Nominal Height	34.0" (86 cm)
Height of Cabinet	33.5" (85 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.63" (60 cm)
Depth with Handle	25.13" (64 cm)
Depth with door at 90°	46.13" (117 cm)
Capacity	3.2 cu.ft. (91 L)
Defrost Type	Manual
Door	White
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Amps	1.3
Weight	105.0 lbs. (48 kg)
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101043500
Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Height Adjustable	1
Interior Height	24.25" (62 cm)
Interior Width	17.75" (45 cm)
Interior Depth	16.5" (42 cm)
Thermostat Type	Digital
Fan Type	Exterior
Refrigerant Type	R134a
Refrigerant Amount	2.7 oz.
High Side PSI	235.0
Low Side PSI	70.0
Level Legs	4
Interior Drawers	Yes
Drawers Qty	3
Compressor Step Height	6.0" (15 cm)
Compressor Step Width	17.75" (45 cm)

Compressor Step Depth	6.0" (15 cm)
Temperature Range	-12°C to -25°C

