

# BAILEY

Model 4124 to 4127

## Economy Electric Hi-Lo Upholstered Mat Tables



**Model 4126**  
Economy Electric  
5 x 7 Mat Table  
[shown with Sand Vinyl  
and Natural Wood Finish]



**Rocker Switch**  
(standard)

**Foot Switch**  
(optional)



### *Model 4124 to 4127 Series Features:*

- 19" to 27" height adjust, using Rocker Switch (standard) or Footswitch (optional).
- Adjustable floor levelers with non-skid rubber feet.
- Herculite Black standard upholstery, with over 15 Naugahyde vinyl colors available at no extra charge.
- Natural Wood Finish Standard (shown above); Optional Dark Walnut Finish available.
- All Hardwood Legs & Frame; clear patient lift clearance on all sides
- Weight capacity 500 lb.
- Ship weight 220 lb. to 330 lb., depending on size.
- Some Assembly Required.

MODEL	SIZE
4124 Upholstered Top	3' x 7'
4125 Upholstered Top	4' x 7'
4126 Upholstered Top	5' x 7'
4127 Upholstered Top	6' x 8'