

BAILEY

Model 3100 Balance Platform

Balance Platform Features:

- 36" x 36" Carpeted Platform.
- Platform height 2" including carpet.
- 500 LB capacity, with solid wood platform, and steel posts and handrails.
- Dual handrails adjust from 28" to 44" high.
- 1 1/2" Diameter Handrails, with Plunger pull and release locking height adjustment
- 30" between uprights allows manual and some power wheelchair access.
- Ships "knocked down", assembly required
- Ships via UPS, weight 89 lbs.

Model 3100
36" x 36" Platform
Covered with HD Low Pile Carpet



Model 3100 Shown With Model 725
Multiple Configurations Possible
Use as mini Training Steps

