

BAILEY

Model 4301 Cabinet Treatment Table

Model 4301 Treatment Table Features:

- 30" height, 30" wide and 78" long; other sizes available.
- 1" thick pressed vinyl upholstery in Camel standard, with your choice of optional Bailey upholstery colors, and either Natural or Dark Walnut wood finishes.
- Full width interior adjustable Shelf.
- Four full access doors in Classic Style fronts; optional Locks available.
- Weight capacity 400 LB
- Ships fully assembled; ship weight 260 LB (Base Model)
- This table is one of our most popular for customized options: Hidden Hinges, Fingerpull Door Edges, add 2" Thick Top, or change to 3 Doors & 3 Drawers; just ask for a quote!



Model 4301
(Shown in Natural, with
Mauve Upholstery)



**Standard
Wire Fingerpulls**



**Model 4301 Customized
No Visible Hardware Option**
(Shown with Custom Wood Finish & 2"
thick Mandarin Orange Upholstery)