

BAILEY

Model 4305 Cabinet Treatment Table

Model 4305 Treatment Table Features:

- 30" height, 30" wide and 78" long.
- 1" thick vinyl upholstery in US 387 Camel standard, with your choice of optional Bailey upholstery colors, and either Natural, Dark Walnut or custom Sherwin Williams wood finishes.
- Six full extension access drawers and One Door in Classic Style fronts; optional Locks available.
- Adjustable Backrest (may be deleted as optional Flat Top)
- Optional Locks, and Professional Style Doors.
- Wire Door and Drawer pulls optional.
- Weight capacity 400 lb.
- Ship weight 275 lb.
- Ships fully Assembled

Model 4305
(Top in Space Blue with Natural
Wood Finish)



Model 4305
(Top in American Beauty with Dark Walnut
Wood Finish)