

# BAILEY

## Models 4310 & 4311 Full Cabinet Style Sports Medicine Treatment Tables

### *Model 4310 & 4311 Treatment Table Common Features:*

- 30" height, 30" wide and 78" long
- 2" thick sewn vinyl upholstery in US 387 Camel standard, with your choice of Bailey upholstery colors, and either Natural or Dark Walnut wood finishes (no cost), or optional Sherwin Williams colors (upcharge).
- Twin Sliding Doors on hardwood tracks with enclosed storage area (Model 4310)
- Four Doors with Wire Fingerpulls (Classic) Style or Optional Recessed Pulls on Pro Style Doors (Model 4311)
- Heavy Duty Adjustable Legrests and Backrest, each with 8 Locking Positions (up to 65 degrees for Backrest, and 32 Degrees for Legrests)
- Weight capacity 350 lb.
- Ships fully assembled; ship weight 260 lb.
- ***Legrests may be ordered as LH or RH; see below examples***



**Model 4311 Cabinet Table**  
Features 4 Door Cabinet,  
Dual Legrests (shown RH) and HD Backrest  
2" Sewn Top shown in American Beauty  
Upholstery, in Natural Wood Finish)



**Model 4310 Cabinet Table**  
Sliding Door Cabinet, Dual Legrests  
(shown LH) and HD Backrest  
2" Sewn Top shown in Imperial Blue  
Upholstery, in Natural Wood Finish)