

BAILEY

Model 4520 Bariatric 4 x 7 Mat Table

Bariatric Mat Table Features:

- 4' W x 7' L x 18"H, with a 2" thick padded top; default vinyl upholstery in US 387 Camel color.
- Choice of over 15 standard Bailey upholstery colors, and Natural Wood finish [other finishes optional].
- ***Six legs and special frame design allows 1000 lb patient capacity!***
- Ships "knocked down", assembly required
- Ship weight 200 lbs.



Bariatric Mat Table
[Shown in US 387 Camel]