

BAILEY

Model 4550 Series BariMatic Electric Hi-Lo Upholstered Mat Tables



***BariMatic Electric Mat Table
[shown in Royal Blue Soft Vinyl]***



***Rocker Switch
(standard)***

Adjustable Feet

Barimatic Mat Table Features:

- 20" to 28" height adjust, using Rocker Switch (standard) (NEW 2024).
- Adjustable floor levelers.
- Camel Soft Vinyl standard upholstery, with optional colors available at no extra charge.
- All steel, heavy duty powder coated frame with crossbrace for stability.
- Lifting Weight capacity 1000 lb.
- Ship weight 270 lb. to 370 lb., depending on size.
- Ships Fully Assembled..

| MODEL | SIZE |
|-------------------------------------|---------|
| 4554 BariMatic with Upholstered Top | 4' x 7' |
| 4555 Barimatic with Upholstered Top | 5' x 7' |
| 4556 Barimatic with Upholstered Top | 6' x 8' |