

BAILEY

Model 497 Couch Recovery Table

Model 497 Features:

- 24" Wide x 20" High x 72" Long, with 2" High Density urethane foam padding.
- Premium knit backed vinyl cover resists wear. Choice of over 15 standard colors [Camel brown is default unless specified].
- Natural Wood Finish Legs Standard, Dark Walnut Finish available. NO plastic laminates.
- Solid Ash Legs, with Solid core Hardwood Plywood construction.
- **Weight Capacity 300 lbs; Ship weight 90 lbs., Ships "knocked down", Assembly required.**
- **Options include Models 46 Paper Holder and 47 Paper Cutter [Assy. Required].**



Model 497
*Shown in Natural Wood Finish
with Peach Upholstery*