


BAILEY

Model 510 Series Floor Mounted Parallel Bars

Model 510 Series Features:

- Height range adjusts from 31" to 41" high, using spring loaded plunger pins.
- Width determined by floor mount spacing; width range is +/- 8 inches.
- Available in 10' thru 24' lengths
- Optional lower Child Handrail available as Model 94A or 94F on all 10 FT & 12 FT bars ONLY.
- Chrome plated and welded carbon steel construction
- **Assembly required**
- **Floor Mounting hardware not provided, to be furnished by installer**
- **5 Year Limited Warranty; weight capacity 300 LB**



Model 510
10 FT Parallel Bars
Floor Mounted
[Height & width Adjustable]

Model 94F
Fixed Height
Child HR Option

| MODEL | LENGTH | POSTS | MODEL | LENGTH | POSTS |
|-------|--------|-------|-------|--------|-------|
| 510 | 10 FT | 4 | 518 | 18 FT | 6 |
| 512 | 12 FT | 4 | 520 | 20 FT | 6 |
| 514 | 14 FT | 6 | 522 | 22 FT | 8 |
| 516 | 16 FT | 6 | 524 | 24 FT | 8 |
| | | | | | |