

BAILEY

Model 530 Series Platform Mounted Parallel Bars, Fixed Width

Features, Model 530 (10' Long) and 532 (12' Long) Series Bars

- Height adjusts from 26" - 44"
- Handrail Width fixed at 25 1/2".
- Width between base posts is 24" for wheelchair access
- Platform size is 30" wide x 120" (530) or 144" (532)
- Chrome Plated Welded Steel Construction with Hardwood Plywood platform
- Positive Plunger Lock Height Adjustment
- Optional Model 94F Child Handrail available
- Optional Model 90 Divider Board available.
- Weight capacity 300 LBS.
- Assembly Required

