

BAILEY

Model 540 Series Platform Mounted Parallel Bars

Features, Models 540-7 (7' Long) 540 (10' Long) and 542 (12' Long) Series Bars

- Height adjusts from 26" - 44"
- Width adjusts from 19" to 25" between bars
- Width between base posts is 36" for wheelchair access
- Platform size is 42" wide x 84" (540-7) 120" (540) or 144" (542)
- Chrome Plated Welded Steel Construction with Hardwood Plywood platform
- Positive Plunger Lock Height Adjustment
- Optional Model 94A or 94F Child Handrail available
- Optional Model 90 Divider Board available (Optional Model 91 Hole Plug is available if no Divider is used)
- Weight capacity 300 LBS.
- Assembly Required

