

BAILEY

Model 595 Wall Mounted Folding Parallel Bars

Model 595 Features:

- Adult Size Parallel Bars that require Wall and Floor Mounting
- Height adjusts from 28" to 42" high
- Width is fixed at 22.5" between handrails and posts
- Chrome plated steel construction; with 1.5" diameter x 84" long handrails with rubber end caps.
- Folds up for storage.
- Weight capacity 250 LB
- Assembly and Wall/Floor Mounting required (mounting hardware not included, provided by installer)



Model 595
Wall Mounted Folding Parallel Bars
[Shown down for operation]



Model 595
Wall Mounted Folding Parallel Bars
[Shown up for storage]