

BAILEY

Model 596 Toddler Folding Parallel Bars

Model 596 Features:

- Suitable for preschool age children.
- Height adjusts 15"/38cm to 21"/53cm, using push lock snap buttons.
- Width adjusts 14"/35.5cm to 23"/58.5cm, by moving base plate bolts.
- Chrome plated steel construction; 1"/25mm diameter x 60"/152cm long handrails with rubber end caps.
- Folds easily for storage.
- Weight capacity 40 lb/18 kg.
- Ship weight 34 lb/15.5 kg.



**Model 596
Toddler Parallel Bars**
*[Shown adjusted to maximum
height and width]*

See our full catalog of products at
www.baileymfg.com