


# BAILEY

## Models 597 (Adult) & 598 (Child) Series Folding Parallel Bars



*Model 597W/598W  
with Wood Bases  
[shown in OPEN position]*



*Model 597S/598S  
with Steel Bases  
[shown in FOLDED position]*

### *Model 597 & 598 Series Features:*

- 31.75" between posts for wheelchair access
- Width between handrails adjusts from 16" to 24"
- 597 Adult Handrail is 1.5" diameter; 598 Child Handrail is 1.0" diameter
- Height adjusts from 22" to 36" (597 Adult) and 19" to 27" (598 Child)
- Base width 42" (Wood and Steel); Wood Base is 1.25" thick, Steel Base is 1/4" thick
- Weight capacity 250 lb. Assembly Required. Folds to 14" Wide and 36" longer folded
- **NOTE: Wood Base models have ONE Stabilizer Bar, and Steel Base Models have TWO Stabilizer Bars.**
- **NOTE: Folding Bars will have noticeable movement while in use. This is normal, and is a design trade off with the ability to fold.**

MODEL	NOTES	MODEL	NOTES
597W-7 Adult, 7 FT L	Wood Bases	597S-7 Adult, 7 FT L	Steel Bases
597W-10 Adult, 7 FT L	Wood Bases	597S-10 Adult, 10 FT L	Steel Bases
598W-7 Child, 7 FT L	Wood Bases	598S-7 Child, 7 FT L	Steel Bases