

BAILEY

Model 6040 Weight Sled

For Occupational Therapy, Work Simulation, Industrial Rehab

Model 6040 Features:

- Versatile work hardening and evaluation tool with four removable weight holders for standard disk weights
- Handle height adjusts from 32" to 44" high in five locking positions
- Plastic Anti-Friction Runners protect floor surfaces. May still damage some floors; provide a non-critical surface (or sheet of plywood) for using this product.
- Accessory Weight Box (Model 6041) available
- Weight capacity 250 lbs/113 Kg.
- Ships via UPS; some assembly required


**Model 6040
Weight Sled**



**Model 6041
Accessory Box**



**Model 6040
uses standard 1" center
hole disk weights**
[shown used on smooth finished
concrete floor]



**Model 6040
with Model 6041
Accessory Box**
[shown used on commercial
unpadded carpet]