

BAILEY

Model 700 Mobile Posture Mirror (Adult)

Model 700 Mobile Mirror Features:

- Adult sized, distortion-free 1/4 inch plate glass mirror
- Glass with Safety Backing per ANSI Z97.1 and 16 CFR 1201
- Full Wood Panel measures 28" x 72" each panel [mirror size 23" x 67"].
- Solid hardwood frame with a durable clear natural finish.
- Heavy duty 2" ball bearing casters feature hard rubber wheels.
- **NOTE: Casters are 2" non-marking Rubber Wheel casters, and may be incompatible with some modern laminate flooring. We are not responsible for marked floors.**



Model 700 Mobile Mirror *[Shown in Natural Wood Finish]*

Total Height 74" / 188cm

*Nominal base Size: 31" / 79cm
wide x 18" / 46cm Deep*

