

BAILEY

Model 701 Mobile Triple Posture Mirror (Adult)

Model 701 Mobile Triple Mirror Features:

- Three Adult sized, distortion-free 1/4 inch plate glass mirrors.
- Glass has Safety Backing per ANSI Z97.1 and 16 CFR 1201.
- Full Hardwood Frame Panels each measure 28" x 72" per panel [mirror size 23" x 67"].
- Solid hardwood frames with a durable clear natural finish.
- Tempered Hardboard back panels.
- **CAUTION:** .Folds into a triangle for moving, and locks into an arc for actual use; do not roll over thresholds or large floor defects.
- **NOTE:** Casters are 2" non-marking Rubber Wheel casters, and may be incompatible with some modern laminate flooring. We are not responsible for marked floors.



**Model 701
Triple Mirror**
*[Shown in
Natural Wood Finish]*

*Total Open Width: 84"
Total Overall Height: 75"
Total Weight: 195 LB*

*Glass Size:
23" / 58cm Wide x
67" / 170cm Tall*

