

BAILEY

Model 702 Wall Posture Mirror (Adult)

Model 702 Mobile Mirror Features:

- Adult sized, distortion-free 1/4 inch plate glass mirror.
- Glass with Safety Backing per ANSI Z97.1 and 16 CFR 1201.
- Full Wood Panel measures 28" x 72" each panel [mirror size 23" x 67"].
- Solid hardwood frame with a durable clear natural finish.
- Tempered Hardboard back panel.
- **NOTE:** Mounting hardware provided by installer to comply with local building codes.



Model 702 Wall Mirror *[Shown in Natural Wood Finish]*

Total Width 28" / 71cm
Total Height 72" / 183cm

Glass Size:
23" / 58cm Wide x
67" / 170cm Tall

