

BAILEY

Model 715 Curb & Ramp Training Set *Wheelchairs - Rehab - Fitness*

Features:

- Four (4) Piece Set of Platform Stools and Ramp
- Non-Slip Alternating Treads with Hardwood Birch Construction & Clear Finish.
- Platform Stool / Curb Sizes: 2" x 32" Square, 4" x 36" square and 6" x 40" square, all with easy pull handles.
- Ramp size 44" wide x 40" long; platform size 44" square, 8" high; overall 84" long.
- Stool / Curbs Nest together for storage under the Ramp.
- Ideal for Wheelchair Training, Step-Ups or beginning Plyometrics
- Weight capacity 300 lb.
- Ships assembled and ready to use.
- *NOTE: The Ramp slope is a non-ADA ramp slope to fit into the allotted size. An ADA compliant version would need to be an overall size of 12 feet long.*

Model 715 Curb & Ramp Set
[Shown with Stools partially unnected for demo]

