

BAILEY

Model 750 & 752 PT Stools

750 & 752 Stool Features:

- 2" Foam Padded 14" x 14" Square Seats
- Chrome Plated Steel Frames
- Both Models are 17" High.
- Seats available in over 15 standard vinyl colors (Camel Brown is default)
- Weight capacity 250 lb.
- Ship weight 14 lb.
- Ships via UPS/Fedex
- Heavy duty 3" ball bearing casters (on 750) feature hard rubber wheels.
- *NOTE: Model 750 Casters are 3" non-marking Rubber Wheel casters, and may be incompatible with some modern laminate flooring. We are not responsible for marked floors.*

