


BAILEY

Model 805, 806 & 807 Convertible Training Stairs

Convertible Staircase Features:

- 24" (805) 30" (806) or 36" (807) wide access, with matching width square Platform.
- 500 LB capacity, with solid hardwood and engineered hardwood construction.
- Triple handrail height design at 18", 27" and 35" for adults and children.
- Powder Coated Steel Upright Handrail Posts
- Step height 6" (4 each) and 3" (8 each), tread depth 10", with anti-slip tread and platform surfaces; 24" high at platform.
- Max height at upper platform handrail is 54".
- Minimum required ceiling height is 101".
- Ships "knocked down", assembly required
- Ship weight 600 lbs. max



Convertible Training Stairs
Straight Configuration



Convertible Training Stairs
Angle Configuration