

# BAILEY

## Model 808 Three Sided Training Stairs

### Staircase Features:

- 30" wide access, with 30" wide x 24" deep Platform at 24"H.
- 500 LB capacity, with solid hardwood and engineered hardwood construction.
- Dual handrail design at 40" and 26" heights for adults and children.
- Step height 6" (4 each), tread depth 10", with anti-slip tread and platform surfaces
- Includes removable Bus Step for public transportation simulation
- Overall Footprint is 55"L x 33"W
- Max height at upper platform handrail is 54".
- Minimum required ceiling height 102" to 108".
- Ships "knocked down", assembly required
- Ship weight 206 lbs.



**808 Training Stairs**  
30" Wide - 3 Sided - "Dead  
End" Configuration



**Bus Step Shown Installed**