

BAILEY

Model 810 Training Stairs 36" Wide - Pass Thru

Staircase Features:

- 36" wide access, with 36" wide x 36" deep Platform at 24"H.
- 500 LB capacity, with solid hardwood and engineered hardwood construction.
- Dual handrail design at 35" and 21" heights for adults and children.
- Step height 6" (8 each), tread depth 10", with anti-slip tread and platform surfaces
- Includes removable Bus Step for public transportation simulation
- Overall Footprint is 96"L x 38"W
- Max height at upper platform handrail is 59".
- Minimum required ceiling height 102" to 108".
- Ships "knocked down", assembly required
- Ship weight 300 lbs.



810 Training Stairs
36" Wide - Pass Thru
Straight Configuration

Bus Step Shown
NOT Installed



Bus Step Shown
Installed