

BAILEY

Model 9505 Electric Tilt Table NEW Revised specifications for 2021

Tilt Table Features:

- Top size is 28" x 78" x 1" thick / 71 cm x 198 cm x 2.5 cm thick
- Height is 31"/78.75 cm
- Walk Off Footplate size is 14" x 17.5" / 35.5 cm x 44.5 cm
- 4 Wheel Locking Swivel Casters, size is 4" / 10 cm diameter
- 1" /2.5 cm thick padded upholstery with Heavyweight vinyl cover. Black is standard color, over 15 other colors are available.
- Three adjustable Black Nylon stabilization/restraint patient straps, with Side Lock Snap Buckles
- Tilts 98°, from vertical to 8° below horizontal
- All Steel Frame is Powder Coated Silver for durability
- NEW - Optional Battery Backup Pack allows operation during power outages
- Weight capacity 350 lb./159 kg.
- Ships fully assembled; ship weight 280 lb.

Model 9505 Electric Tilt Table

NEW Drive System design
NEW low voltage Handswitch
NEW higher patient capacity
NEW Optional Battery Pack

