

Bariatric Sentra EC Heavy-Duty, Extra-Extra-Wide Wheelchair

HCPCS: K0007-1,K0007-1



- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable reinforced nylon upholstery
- Steel, spoke wheels are strong and maintenance free
- Heavy-duty, flat-free 8" front casters are adjustable in two positions
- Extra-heavy-duty front forks
- Reinforced side frame, caster journals and dual cross bars
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Extra-large, padded armrests provide comfort and support
- Riggings have extra-large, heavy-duty, aluminum foot plates that provide maximum support, are attractive and crack proof (Must be ordered separately, PH-SF for swing-away footrests, PHELR for elevating legrests)
- Comes with push-to-lock wheel locks
- Weight Capacity 700 lbs

PRODUCTS	Description
STD26ECDDA	Sentra EC,26"Det.DeskArm,1/cs
STD26ECDFA	Sentra EC,26in Det.FullArm,1/cs
STD28ECDDA	Sentra EC,28"Det.DeskArm,1/cs
STD28ECDFA	Sentra EC,28"Det.FullArm,1/cs
STD30ECDDA	Sentra EC,30"Det.DeskArm,1/cs
STD30ECDFA	Sentra EC,30"Det.FullArm,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	62 lb 64 lb 66 lb
Product Dimensions	43" (L) w/ riggings x 34.5" (W) x 35" (H)
Seat Dimensions	26" (W) x 20" (D)
Weight Capacity	700 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	68 lb 70 lb 72 lb
Folded Dimensions	43" (L) w/ riggings x 14" (W) x 35" (H)
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov]