



CrossRunner™ Soft Knee Brace



The CrossRunner is a soft-hinged knee brace designed to provide optimal support with maximum comfort in every detail. The perfect cross between comfort and support.

Product Features

- Universal fit to accommodate both left and right knees
- Composed of a hybrid fabric, BreathFit™, that allows for both breathability and compression
- Adjustable range of motion (ROM) control with snap-fit ROM stops
- Durable design of both the sleeve and hinge
- Wraparound version features flaps for a secure fit
- Pull-on version features pull-tabs for easy application

Common Examples of Use

- Most ACL, PCL, MCL, and LCL injuries and instabilities
- Meniscal instabilities
- Patellar instabilities

- Mild osteoarthritis

Categories: Knee Braces, Adjustable ROM, Soft Supports.

Instructions for Use

[Wraparound Instructions for Use](#)

[Pull-On Instructions for Use](#)

Part Numbers & Sizing Chart

Wraparound

Standard Part #	Short Part #	Size	Thigh Circumference	Calf Circumference
KN1000-005	KN3000-005	XXS	12" - 14" (30.5 - 35.5 cm)	9.5" - 10.5" (24.1 - 26.7 cm)
KN1000-010	KN3000-010	XS	14" - 16" (35.5 - 40.6 cm)	10.5" - 11.5" (26.7 - 29.2 cm)
KN1000-020	KN3000-020	S	16" - 18" (40.6 - 45.7 cm)	11.5" - 12.5" (29.2 - 31.8 cm)
KN1000-030	KN3000-030	M	18" - 21" (45.7 - 53.3 cm)	12.5" - 15" (31.8 - 38.1 cm)
KN1000-040	KN3000-040	L	21" - 24" (53.3 - 60.9 cm)	15" - 18" (38.1 - 45.7 cm)
KN1000-050	KN3000-050	XL	24" - 27" (60.9 - 68.5 cm)	18" - 20" (45.7 - 50.8 cm)
KN1000-060	KN3000-060	2XL	27" - 31" (68.5 - 78.7 cm)	20" - 23.5" (50.8 - 59.7 cm)
KN1000-070	KN3000-070	3XL	31" - 36" (78.7 - 91.4 cm)	23.5" - 27" (59.7 - 68.6 cm)

Pull-On

Standard Part #	Short Part #	Size	Thigh Circumference	Calf Circumference
KN2000-010	KN4000-010	XS	14" - 16" (35.5 - 40.6 cm)	10.5" - 11.5" (26.7 - 29.2 cm)
KN2000-020	KN4000-020	S	16" - 18" (40.6 - 45.7 cm)	11.5" - 12.5" (29.2 - 31.8 cm)
KN2000-030	KN4000-030	M	18" - 21" (45.7 - 53.3 cm)	12.5" - 15" (31.8 - 38.1 cm)
KN2000-040	KN4000-040	L	21" - 24" (53.3 - 60.9 cm)	15" - 18" (38.1 - 45.7 cm)
KN2000-050	KN4000-050	XL	24" - 27" (60.9 - 68.5 cm)	18" - 20" (45.7 - 50.8 cm)
KN2000-060	KN4000-060	2XL	27" - 31" (68.5 - 78.7 cm)	20" - 23.5" (50.8 - 59.7 cm)
KN2000-070	KN4000-070	3XL	31" - 36" (78.7 - 91.4 cm)	23.5" - 27" (59.7 - 68.6 cm)

Standard length is recommended for a person that is 5' 5" (165 cm) and over.

Thigh circumference measured at 6" (15 cm) above center of patella.

Calf circumference measured at 6" (15 cm) below center of patella.

Part #	Description
KN0000-001	Kit, Condyle and Hook, Pair
KN0000-002	Kit, Patella Buttress