

# ACR82L



61.25" x 21" x 21.88" (H x W x D)

7.56 cu.ft. general purpose pharmaceutical all-refrigerator with self-closing glass door and digital controls

**Highlights:**

Made in North America in a TAA designated country

Adjustable cantilevered shelving for flexible storage

Self-closing door ensures stored items are not accidentally overexposed to warm temperatures


**Product Features:**

General purpose healthcare refrigerator made in North America	All-refrigerator designed to operate between 2 and 8°C for a variety of storage applications including cold packs, laboratory specimens, and medications
Temperature Range: 2° and 8°C	Preset to 5°C as recommended for most items stored in general purpose applications
TAA compliant	Ideal for GSA contracts and other applications requiring a unit to have been manufactured in a TAA designated country
Glass door	Tempered glass door provides heat-safe view of interior contents
Self-closing door	Door automatically closes if left ajar to help protect items in high traffic settings
Pocket handle	Handle recessed in the side of the door for more efficiency in crowded facility settings
Fan-cooled interior	Powerful fan-forced cooling to ensure fast temperature recovery
Automatic defrost	Reduced user maintenance with automatic defrost
High/low temperature alarm	Audible and visual alarm warns if the temperature goes out of range


ACR82L Specifications:

Overview	
Nominal Height	61.0" (155 cm)
Height of Cabinet	60.88" (155 cm)
Height to Hinge Cap	61.25" (156 cm)
Nominal Width	21.0" (53 cm)
Width	21.0" (53 cm)
Nominal Depth	22.0" (56 cm)
Depth	21.88" (56 cm)
Depth with door at 90°	41.5" (105 cm)
Capacity	7.56 cu.ft. (214 L)
Defrost Type	Automatic
Door	Glass
Cabinet	White
US Electrical Safety	UL
Canadian Electrical Safety	ULC
Sanitation	UL-S
Amps	2.5
Voltage/Frequency	115 V AC/60 Hz
Weight	170.0 lbs. (77 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101086071
Refrigerator Features	
Door Swing	RHD
Reversible	No
Shelf Type	Plastic
Shelf Qty	3
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R290
Refrigerant Amount	2.1oz.
High Side PSI	293.0
Low Side PSI	135.0
Interior Light	Yes
Temperature Range	2° to 8°C
Sound Level	44.3 dB
Dimensions	

Interior Height	39.38" (100 cm)
Interior Width	17.5" (44 cm)
Interior Depth	18.25" (46 cm)



TOP VIEW



FRONT VIEW