

ADA305AFSSHVLHD



32" x 20.88" x 22.5" (H x W x D)

Built-in undercounter -25°C ADA compliant all-freezer with stainless steel front and handle, factory installed lock, digital controls, interior baskets, hospital cord with 'green dot' plug, access port, and manual defrost operation

Highlights:

Flexible design allows built-in or freestanding use in 21" wide spaces

Easy temperature management with a digital thermostat

Basket shelving offers convenient storage of a variety of items

Product Features:

20" Wide Built-In Healthcare All-Freezer, ADA-Compliant	Slim-fitting all-freezer for general healthcare applications at an ADA-compliant height
ADA-Compliant Design	32" height allows easy installation under lower ADA compliant counters
-25°C Capable	Low temperature operation ideal for vaccines, pharmaceuticals, and other sensitive materials
Built-In Capable	Front-breathing design lets you make the best use of space by installing your appliance under the counter
Stainless Steel Wrapped Door	Durable and stylish design with lasting elegance
Professional Stainless Steel Handle	Thin (14mm diameter) stainless steel handle mounted vertically for a sleek look that matches many popular styles of cabinetry
Factory-Installed Lock	Located towards the top of the door for easy access and to provide added security of unit (2 keys included)
Access Port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Adjustable Digital Thermostat	Electronic controls located in the front kickplate for easy & intuitive temperature management, with a readout in Fahrenheit or Celsius

ADA305AFSSHVLHD Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	31.75" (81 cm)
Height to Hinge Cap	32.0" (81 cm)
Nominal Width	21.0" (53 cm)
Width	20.88" (53 cm)
Nominal Depth	23.0" (58 cm)
Depth	22.5" (57 cm)
Depth with Handle	24.0" (61 cm)
Depth with door at 90°	40.25" (102 cm)
Capacity	2.47 cu.ft. (70 L)
Defrost Type	Manual
Door	Stainless Steel
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	75.0 lbs. (34 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101151328
Freezer	
Door Swing	LHD
Reversible	Yes
Interior Height	23.75" (60 cm)
Interior Width	14.5" (37 cm)
Interior Depth	17.0" (43 cm)
Shelf Type	Drawers
Shelf Qty	3
Thermostat Type	Digital
Fan Type	Exterior
Refrigerant Type	R600a
Refrigerant Amount	1.34 oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs	2
Interior Drawers	Yes
Drawers Qty	3

Compressor Step Height	5.75" (15 cm)
Compressor Step Width	14.5" (37 cm)
Compressor Step Depth	6.0" (15 cm)
Temperature Range	-25° to -15°C

