

ADA305AFSSTBC



32" x 19.75" x 22.5" (H x W x D)

Built-in undercounter -25°C ADA compliant all-freezer in white with a stainless steel door, lock, pure copper towel bar handle, digital controls, interior baskets, hospital cord with 'green dot' plug, factory installed access port, and manual defrost operation

Highlights:

Medical grade design to function at lower temperatures (-25°C) than standard residential units

Just 32" high for an easy fit under lower ADA compliant counters

Copper towel bar handle kills microbes in a shorter timeframe than other metals for added safety

Product Features:

ADA compliant design	32" height allows easy installation under lower ADA compliant counters
-25°C capable	Low temperature operation ideal for vaccines, pharmaceuticals, and other sensitive materials
Built-in capable	Front-breathing design lets you make the best use of space by installing your appliance under the counter
Copper handle	Naturally antimicrobial, the towel bar handle made from pure copper helps to protect against microbes for added safety and cleanliness
Factory-installed lock	Located towards the top of the door for easy access and to provide added security of unit (2 keys included)
Access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Adjustable digital thermostat	Electronic controls located in the front kickplate for easy & institutive temperature management, with a readout in Fahrenheit or Celsius
Basket shelves	Three slide-out baskets offer convenient storage, with plastic fronts to help protect contents from warmer air when the door is opened

3" thick door Improved insulation for an efficient performance; NOTE: due to the oversized design, the door does not open in its own radius. If installing this freezer against a wall, please see our drawings to ensure adequate clearance. SUMMIT APPLIANCE DIVISION, FELIX STORCH, INC. ISO 9001:2015 Certified 770 Garrison Avenue Bronx, NY 10474 USA TEL 718-895-1100 Einfo@summitappliance.com  www.summitappliance.com

ADA305AFSSTBC Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	31.75" (81 cm)
Height to Hinge Cap	32.0" (81 cm)
Nominal Width	20.0" (51 cm)
Width	19.75" (50 cm)
Width with Door Open	20.25" (51 cm)
Nominal Depth	23.0" (58 cm)
Depth	22.5" (57 cm)
Depth with Handle	25.0" (64 cm)
Depth with door at 90°	40.25" (102 cm)
Capacity	2.47 cu.ft. (70 L)
Defrost Type	Manual
Door	Stainless Steel
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	75.0 lbs. (34 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101077734
Freezer	
Door Swing	RHD
Reversible	Yes
Interior Height	23.75" (60 cm)
Interior Width	14.5" (37 cm)
Interior Depth	17.0" (43 cm)
Shelf Type	Drawers
Shelf Qty	3
Thermostat Type	Digital
Fan Type	Exterior
Refrigerant Type	R600a
Refrigerant Amount	1.34 oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs	2
Interior Drawers	Yes

Drawers Qty	3
Compressor Step Height	5.75" (15 cm)
Compressor Step Width	14.5" (37 cm)
Compressor Step Depth	6.0" (15 cm)
Temperature Range	-25° to -15°C

