

ADA404REF



32" x 19.75" x 22.5" (H x W x D)

Built-in undercounter ADA compliant +2°C to +8°C all-refrigerator in white with lock, digital controls, wire shelving, hospital cord with 'green dot' plug, factory installed access port, and automatic defrost operation

Highlights:

Just 32" high for an easy fit under lower ADA compliant counters

Flexible design allows built-in or freestanding use in 20" wide spaces

Wire shelving offers convenient storage of a variety of items

Product Features:

ADA compliant design	32" height allows easy installation under lower ADA compliant counters
Built-in capable	Front-breathing design lets you make the best use of space by installing your appliance under the counter
2 to 8° operation	Designed to safely store select pharmaceuticals and other materials with strict temperature requirements
Factory-installed lock	Located towards the top of the door for easy access and to provide added security of unit (2 keys included)
Access port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Adjustable digital thermostat	Electronic controls located in the front kickplate for easy & institutive temperature management, with a readout in Fahrenheit or Celsius
Adjustable wire shelves	Four wire shelves provide optimal storage and improved air circulation in the 3.32 cu. ft. interior
3" thick door	Improved insulation for an efficient performance; NOTE: due to the oversized design, the door does not open in its own radius. If installing this freezer against a wall, please see our drawings to ensure adequate space

ADA404REF Specifications:

Overview

Nominal Height	32.0" (81 cm)
Height of Cabinet	31.75" (81 cm)
Height to Hinge Cap	32.0" (81 cm)
Nominal Width	20.0" (51 cm)
Width	19.75" (50 cm)
Width with Door Open	20.25" (51 cm)
Nominal Depth	23.0" (58 cm)
Depth	22.5" (57 cm)
Depth with Handle	24.38" (62 cm)
Depth with door at 90°	40.25" (102 cm)
Capacity	3.32 cu.ft. (94 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	59.0 lbs. (27 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101063058

Refrigerator Features

Door Swing	RHD
Reversible	Yes
Shelf Type	Wire
Shelf Qty	4
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.99oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs Qty	2
Interior Light	Yes
Temperature Range	2 to 8°C

Dimensions

Interior Height	23.5" (60 cm)
-----------------	---------------

Interior Width	16.25" (41 cm)
Interior Depth	17.0" (43 cm)
Compressor Step Height	6.25" (16 cm)
Compressor Step Width	16.25" (41 cm)
Compressor Step Depth	6.0" (15 cm)

