

ADA404REFCAL



32" x 19.75" x 22.5" (H x W x D)

Built-in undercounter ADA compliant +2°C to +8°C all-refrigerator in white with lock, digital controls, wire shelving, hospital cord with 'green dot' plug, NIST calibrated thermometer with alarm, and automatic defrost operation

Highlights:

- Just 32" high for an easy fit under lower ADA compliant counters
- Flexible design allows built-in or freestanding use in 20" wide spaces
- High temperature alarm includes NIST calibrated external temperature readout (NIST certificate included)

Product Features:

ADA compliant design	32" height allows easy installation under lower ADA compliant counters
Built-in capable	Front-breathing design lets you make the best use of space by installing your appliance under the counter
+2°C to +8°C operation	Designed to safely store vaccines, select pharmaceuticals, and other materials with strict temperature requirements
NIST calibrated external temperature display	Alarm includes continuous NIST calibrated display of the current interior and high/low temperature, and can be switched from Celsius to Fahrenheit at the touch of a button
Certificate of calibration included	Calibrated to NIST standards in our ISO/IEC 17025:2017 accredited laboratory, with a certificate included for your records
Buffered temperature probe	Thermometer probe is encased in a glycol-filled bottle and stored in a fixed box on the left wall of the inner cabinet. This is the ideal way to produce temperature readings closest to stored contents, as recommended by the CDC, VFC and VFA's guidelines
3" thick door	Improved insulation for an efficient performance; NOTE: due to the oversized design, the door does not open in its own radius. If installing this freezer against a wall, please see our drawings to ensure adequate space

Flat door liner Interior design conforms to CDC guidelines to ensure stable conditions for sensitive materials

ADA404REFCAL Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	31.75" (81 cm)
Height to Hinge Cap	32.0" (81 cm)
Nominal Width	20.0" (51 cm)
Width	19.75" (50 cm)
Width with Door Open	20.25" (51 cm)
Nominal Depth	23.0" (58 cm)
Depth	22.5" (57 cm)
Depth with Handle	24.38" (62 cm)
Depth with door at 90°	40.25" (102 cm)
Capacity	3.32 cu.ft. (94 L)
Defrost Type	Automatic
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	59.0 lbs. (27 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101079059
Refrigerator Features	
Door Swing	RHD
Reversible	Factory Reversible
Shelf Type	Wire
Shelf Qty	4
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.99oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs Qty	2
Interior Light	Yes
Sound Level	2°C to 8°C
Dimensions	

Interior Height	23.5" (60 cm)
Interior Width	16.25" (41 cm)
Interior Depth	17.0" (43 cm)
Compressor Step Height	6.25" (16 cm)
Compressor Step Width	16.25" (41 cm)
Compressor Step Depth	6.0" (15 cm)

