

ADA404REFIF



32" x 20.25" x 23" (H x W x D)

Built-in undercounter ADA compliant +2°C to +8°C all-refrigerator with with a panel-ready door, , factory-installed lock, digital controls, wire shelving, hospital cord with 'green dot' plug, and automatic defrost operation

Highlights:

Just 32" high for an easy fit under lower ADA compliant counters

Flexible design allows built-in or freestanding use in 21" wide spaces

Wire shelving offers convenient storage of a variety of items

Product Features:

20" Wide Built-In Healthcare All-Refrigerator, Slim-fitting all-refrigerator for general healthcare applications at an ADA-compliant height ADA-Compliant	
Built-In Capable	Front-breathing design lets you make the best use of space by installing your appliance under the counter
ADA-Compliant Design	32" height allows easy installation under lower ADA compliant counters
Panel-Ready Door	Integrated stainless steel door frame accepts full overlay panels, allowing you to customize the refrigerator door or blend it into your cabinetry design; panel size:19 13/16" W x 27 1/8" H; customers provide their own handles
2 to 8°C Operation	Designed to safely store select pharmaceuticals and other materials with strict temperature requirements
Automatic Defrost	Reduced user maintenance with an auto defrost system
Access Port	Unit comes with a factory installed probe hole for users to connect their own monitoring equipment
Factory-Installed Lock	Located towards the top of the door for easy access and to provide added security of unit (2 keys included)

ADA404REFIF Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	31.75" (81 cm)
Height to Hinge Cap	32.0" (81 cm)
Nominal Width	21.0" (53 cm)
Width	20.25" (51 cm)
Nominal Depth	23.0" (58 cm)
Depth	23.0" (58 cm)
Depth with Handle	25.0" (64 cm)
Capacity	3.32 cu.ft. (94 L)
Defrost Type	Automatic
Door	Panel-Ready
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	75.0 lbs. (34 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101151311
Refrigerator Features	
Door Swing	RHD
Reversible	Factory Reversible
Height Adjustable	No
Shelf Type	Wire
Adjustable Shelves	Yes
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	0.99oz.
High Side PSI	180.0
Low Side PSI	80.0
Level Legs Qty	4
Interior Light	Yes
Temperature Range	2 to 8°C
Dimensions	
Interior Height	23.5" (60 cm)

Interior Width	16.25" (41 cm)
Interior Depth	17.0" (43 cm)
Compressor Step Height	6.25" (16 cm)
Compressor Step Width	16.25" (41 cm)
Compressor Step Depth	6.0" (15 cm)

