

ARG31PVBIADA-

KLHD 64" x 10.75" x 22.63" (H x W x D)



PHARMA-VAC Performance Series stacked combination of ARG31PVBIADA vaccine all-refrigerator and AFZ2PVBIADA manual defrost all-freezer

Highlights:

- Purpose-built design to store medication, pharmaceuticals, and other sensitive materials
- Independent refrigerator and freezer storage in one slim-fitting footprint for space-challenged facilities
- Intelligent microprocessor digital temperature controller for optimum temperature management

Product Features:

Performance Series Refrigerator-Freezer Combination	Purpose-built all-refrigerator and all-freezer for vaccine storage in medical, pharmaceutical, and life science applications
Purpose-built design	Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC/VFA guidelines
+/-1°C Variation	The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period
Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
Easy assembly	Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top
Antimicrobial silver ion-coated handles	Powder-coated handles with naturally occurring silver ions to help reduce the spread of germs
Adjustable temperature	Intuitive controls let you adjust the interior temperature between +2 to +8°C (refrigerator) and -30° to -10°C (freezer)
Temperature alarms	Visual and audible alarms alert if the interior temperature rises or falls out of set range

ARG31PVBIADA- AFZ2PVBIADASTACKLHD Specifications:

Overview	
Nominal Height	64.0" (163 cm)
Height of Cabinet	64.0" (163 cm)
Nominal Width	20.0" (51 cm)
Width	19.75" (50 cm)
Width with Door Open	20.25" (51 cm)
Nominal Depth	23.0" (58 cm)
Depth	22.63" (57 cm)
Depth with Handle	24.38" (62 cm)
Depth with door at 90°	41.25" (105 cm)
Capacity	5.3 cu.ft. (150 L)
Door	Glass
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Voltage/Frequency	115 V AC/60 Hz
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101104126
Refrigerator-Freezer	
Door Swing	LHD
Reversible	No
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	24.0" (61 cm)
Refrigerator - Interior Width	16.0" (41 cm)
Refrigerator - Interior Depth	15.0" (38 cm)
Refrigerator - Capacity	2.83 cu.ft
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	4
Refrigerator - Temperature Range	2 to 8°C
Refrigerator - Refrigerant Type	R600a
Refrigerator - Refrigerant Amount	1.41 oz.
Freezer - Interior Height	22.0" (56 cm)
Freezer - Interior Width	14.25" (36 cm)

Freezer - Interior Depth	15.5" (39 cm)
Freezer - Capacity	2.47 cu.ft
Freezer - Shelf Type	Wire
Freezer - Shelf Qty	4
Thermostat Type	Digital
Fan Type	Interior
Freezer - Temperature Range	-30 to -10°C
Freezer - Refrigerant Type	R290a
Freezer - Refrigerant Amount	1.31 oz.

