

ARG6PV-AFZ1PVSTACKLHD



54.75" x 23.38" x 24.38" (H x W x D)

PHARMA-VAC Performance Series stacked combination of ARG6PV vaccine all-refrigerator and AFZ1PV manual defrost all-freezer

Highlights:

- Purpose-built design to store medication, pharmaceuticals, and other sensitive materials
- Independent refrigerator and freezer storage in one slim-fitting footprint for space-challenged facilities
- Intelligent microprocessor digital temperature controller for optimum temperature management

Product Features:

Performance Series Refrigerator-Freezer Combination	Purpose-built all-refrigerator and all-freezer for vaccine storage in medical, pharmaceutical, and life science applications
Purpose-built design	Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC/VFA guidelines
+/-1°C Variation	The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period
Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
Easy assembly	Unit ships in three boxes for customer to install the StackRack on the bottom refrigerator/freezer and place and secure the other unit on top
Antimicrobial silver ion-coated handles	Powder-coated handles with naturally occurring silver ions to help reduce the spread of germs
Adjustable temperature	Intuitive controls let you adjust the interior temperature between +2 to +8°C (refrigerator) and -30° to -10°C (freezer)
Temperature alarms	Visual and audible alarms alert if the interior temperature rises or falls out of set range

ARG6PV-AFZ1PVSTACKLHD Specifications:

Overview	
Nominal Height	55.0" (140 cm)
Height of Cabinet	54.75" (139 cm)
Nominal Width	24.0" (61 cm)
Width	23.38" (59 cm)
Nominal Depth	25.0" (64 cm)
Depth	24.38" (62 cm)
Capacity	7.4 cu.ft. (210 L)
Door	Glass
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Voltage/Frequency	115 V AC/60 Hz
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101145921
Refrigerator-Freezer	
Door Swing	LHD
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	24.5" (62 cm)
Refrigerator - Interior Width	20.0" (51 cm)
Refrigerator - Interior Depth	17.0" (43 cm)
Refrigerator - Capacity	6.0 cu.ft
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	4
Refrigerator - Temperature Range	2 to 8°C
Refrigerator - Refrigerant Type	R600a
Refrigerator - Refrigerant Amount	1.66 oz.
Freezer - Interior Height	16.5" (42 cm)
Freezer - Interior Width	15.25" (39 cm)
Freezer - Interior Depth	15.75" (40 cm)
Freezer - Capacity	1.4 cu.ft
Freezer - Shelf Type	Wire
Freezer - Shelf Qty	2

Thermostat Type	Digital
Fan Type	Interior
Freezer - Temperature Range	-30 to -10°C
Freezer - Refrigerant Type	R290a
Freezer - Refrigerant Amount	0.85 oz.

